


# <sup>13</sup>C-Enriched Carbon Sources for Solid State NMR of Proteins

• BioExpress® 1000 • Glycerol • Glucose



## Tailoring <sup>13</sup>C Enrichment of Proteins for Solid State NMR

For many years, the use of highly enriched, uniformly isotope-labeled recombinant proteins has enabled the study of protein structure and dynamics by NMR spectroscopy. This approach has been relatively easy to perform since uniformly <sup>13</sup>C-labeled carbon sources such as <sup>13</sup>C<sub>6</sub> glucose have been readily available. Unfortunately, uniform <sup>13</sup>C labeling of highly enriched proteins is not conducive to some NMR experiments (such as solid state NMR experiments using stationary aligned samples) because strong, homonuclear <sup>13</sup>C-<sup>13</sup>C dipole-dipole coupling can act to broaden <sup>13</sup>C resonances very significantly, sometimes so severely that they are beyond detection. An effective approach to alleviating this problem is to use proteins containing spatially isolated <sup>13</sup>C sites. Such sparsely labeled proteins can either be produced using uniform, random, fractionally enriched carbon sources or specific <sup>13</sup>C-labeled carbon sources.

## Uniform Fractionally <sup>13</sup>C-Enriched Bacterial Cell Growth Media and Glucose (<sup>13</sup>C<sub>6</sub>)

CIL is also pleased to offer fractionally <sup>13</sup>C-enriched BioExpress® 1000, as well as fractionally enriched glucose, for those researchers that use minimal essential media instead of a ready-made rich culture medium.

Catalog No.	Description	Amount
CGM-1000-CN-25-S	BioExpress® 1000 ( <sup>13</sup> C, 25%; <sup>15</sup> N, 98%)	10 mL (10x)
CGM-1000-CN-25	BioExpress® 1000 ( <sup>13</sup> C, 25%; <sup>15</sup> N, 98%)	100 mL (10x)
CGM-1000-CN-35-S	BioExpress® 1000 ( <sup>13</sup> C, 35%; <sup>15</sup> N, 98%)	10 mL (10x)
CGM-1000-CN-35	BioExpress® 1000 ( <sup>13</sup> C, 35%; <sup>15</sup> N, 98%)	100 mL (10x)
CLM-1396-25-1	Glucose ( <sup>13</sup> C, 24-25%)	1 g

The number of adjacent <sup>13</sup>C sites in most amino acid residues can be minimized by using uniform fractionally <sup>13</sup>C-enriched bacterial cell growth media or glucose [<sup>13</sup>C<sub>6</sub>]. Opella et al.<sup>1</sup> have reported optimal results using BioExpress® media with fractional <sup>13</sup>C enrichments between 25% and 35%.

## Reference

1. Filipp, F.V.; Sinha, N.; Jairam, L.; Bradley, J.; Opella, S.J. **2009**. Labeling strategies for <sup>13</sup>C-detected aligned-sample solid-state NMR of proteins, *J Magn Reson*, 201, 121-130.

## Specific <sup>13</sup>C-Labeled Carbon Sources

CIL offers the following specific <sup>13</sup>C-labeled substrates for use in the microbial expression of proteins labeled with isolated <sup>13</sup>C sites.

Catalog No.	Description	Amount
CLM-1397	Glycerol (2- <sup>13</sup> C, 99%)	1g
CLM-1857	Glycerol (1,3- <sup>13</sup> C, 99%)	1g
CLM-746	Glucose (2- <sup>13</sup> C, 99%)	1g
CLM-441	Sodium bicarbonate ( <sup>13</sup> C, 99%)	1g

Glycerol [1,3-<sup>13</sup>C<sub>2</sub>] and glycerol [2-<sup>13</sup>C] have been shown to give rise to complementary labeling patterns. When choosing the appropriate specific <sup>13</sup>C-labeled carbon source, the <sup>13</sup>C enrichment of the expressed protein can be tailored to highlight regions of interest. Sodium bicarbonate (<sup>13</sup>C) should be added to growths using 2-<sup>13</sup>C glycerol.

## Advantages

### Glycerol (2-<sup>13</sup>C)

- labels alpha carbons
- increases overall labeling of carbons
- significantly lower levels of carbonyl and aliphatic sidechain carbons

### Glycerol (1,3-<sup>13</sup>C<sub>2</sub>)

- creates reduced intensity in the α-carbon region.<sup>1</sup>

### Glucose (2-<sup>13</sup>C)

- minimizes adjacent <sup>13</sup>C pairs in most residues
- reported to be more useful than specifically labeled glycerol in detecting backbone <sup>13</sup>C resonances.<sup>1</sup>