

Stable Isotope-Labeled Mixtures, Sets, and Kits

For Mass Spectrometric Applications




Table of Contents

Mixtures and Sets		Kits
Amino Acid Mixes	3	Metabolo
Bile Acid Mixes	7	Metaboli ⁻
Carnitine/Acylcarnitine Mixes	8	IROA® Bio
Fatty Acid/Lipid Mixes	11	PeptiQua
Nucleic Acid Mixes	12	INLLIGHT
Organic Acid Mixes	14	
Other Mixes	16	Q & As
Steroid Mixes	19	

Kits	
Metabolomic QC Kits	21
Metabolite Extract Kits	24
IROA® Biochemical Quantitation Kits	29
PeptiQuant™ Plus Assay Kits	31
INLLIGHT® Glycan Tagging Kit	34
Q & As	35

Please visit **isotope.com** for pricing and availability.

For research use only. Not for use in diagnostic procedures.

The global and targeted measurement of biomolecules continues to be two areas of growing focus in analytical chemistry and biomedicine. Drivers for this research include efforts to better understand the underlying mechanisms of disease pathogenesis and to improve precision medicine through the qualitative/quantitative screening of candidate biomarkers. To address these objectives, mass spectrometry (MS)-based approaches are increasingly utilized and have been aided by advancements in experimental methodologies, instrumentation technologies, and bioinformatic tools.

To help advance research in MS 'omics and MS/MS screening, Cambridge Isotope Laboratories, Inc. (CIL) is pleased to offer the largest variety of stable isotope-labeled sets, mixtures, and kits. The mixes are formulated neat and/or as solutions and are readily available for immediate use. In addition to the packaged mix(es), the kits include a user manual and shipping documents (e.g., certificate of analysis, safety data sheet). The manual outlines general procedures and processing examples for user reference, as well as troubleshooting notes and analysis results. Please see isotope.com for product details, pricing, and to inquire about customized mixes.

Benefits

- easy to implement
- offers end user flexibility
- reduces development time and cost
- enhances data quality
- improves reproducibility
- renders confidence in analytical results

Features Overview

- mixes supplied in neat and/or solution form
- preferentially ¹³C- and/or ¹⁵N-enriched
- site-specific or uniform labeling
- established specification guidelines
- procedural guides in kit manuals
- broad-spanning applications (from QC to quantification)



Mixtures and Sets

The mixtures available off-the-shelf are class-specific for 'omics and MS/MS screening applications. These mixes are amenable for use in quality control and qualification/quantification exercises using targeted, semi-targeted, or untargeted LC-MS methodologies. Outlined below is an overview of our current mix offerings, as well as details into their compositions and usage specifications (i.e., reconstitution guidelines, storage, and stability). For reference purposes, example results and published manuscripts are also provided.

Amino Acid Mixes and Sets

Amino acids (AAs) play critical roles in biological functions as both building blocks of peptides/proteins and intermediates of various metabolic pathways (e.g., citric acid cycle, urea cycle). These compounds are also reported to influence the pathogenesis and propagation of metabolic disorders/disease. To aid continued development and application, CIL has formulated a number of stable isotope-labeled (and unlabeled) AA mixtures. These include mixes of the canonical amino acids (MSK-A2 and MSK-CAA), the rare or unnatural noncanonical amino acids (MSK-NCAA), and a series of reference standard AA mixes (e.g., NSK-A).

Overview

Catalog No.	Description	No. of Metabolites	Unit Size
MSK-A2	Metabolomics Amino Acid Mix	17	1.2 mL
MSK-CAA	Canonical Amino Acid Mix	20	1 vial
MSK-NCAA	Noncanonical Amino Acid Mix	7	1 vial
NSK-A	Amino Acid Standard Mix Set A	12	1 vial, 10 vials
NSK-A1	Amino Acid Standard Mix Set A1	12	1 vial, 10 vials
NSK-A-TS	Amino Acid Tuning Standard Mix Set A	5	1 vial
NSK-AB	Standard Mix Sets A and B	12 (in A), 8 (in B)	2 × 10 vials





Click on the thumbnails or visit isotope.com/applications/ for more information.

MSK-A2 and MSK-CAA Mixes

Composition

The A2 mix (in 0.1 M HCl) consists of 17 amino acids (in solution), and the CAA mix (dried down) comprises 20. Reconstituting the CAA mix in 1 mL solvent (e.g., 50% methanol) results in 2.5 mM concentrations (exception: L-cystine at 1.25 mM). **Note:** This CAA mix should not be reconstituted in 0.1 M HCl as Asn and Gln are unstable in acid.

Compound	Abbrev.	Label and Enrichment	Conc. (mM)
L-Alanine	Ala	¹³ C ₃ , 99%; ¹⁵ N, 99%	2.5
L-Arginine·HCl	Arg	¹³ C ₆ , 99%; ¹⁵ N ₄ , 99%	2.5
L-Asparagine·H ₂ O*	Asn	¹³ C ₄ , 99%; ¹⁵ N ₂ , 99%	2.5
L-Aspartic acid	Asp	¹³ C ₄ , 99%; ¹⁵ N, 99%	2.5
L-Cystine	Cys-Cys	¹³ C ₆ , 99%; ¹⁵ N ₂ , 99%	1.25
L-Glutamic acid	Glu	¹³ C ₅ , 99%; ¹⁵ N, 99%	2.5
L-Glutamine*	Gln	¹³ C ₅ , 99%; ¹⁵ N ₂ , 99%	2.5
Glycine	Gly	¹³ C ₂ , 99%; ¹⁵ N, 99%	2.5
L-Histidine·HCl·H ₂ O	His	¹³ C ₆ , 97-99%; ¹⁵ N ₃ , 97-99%	2.5
L-Isoleucine	lle	¹³ C ₆ , 99%; ¹⁵ N, 99%	2.5
L-Leucine	Leu	¹³ C ₆ , 99%; ¹⁵ N, 99%	2.5
L-Lysine-2HCl	Lys	¹³ C ₆ , 99%; ¹⁵ N ₂ , 99%	2.5
L-Methionine	Met	¹³ C ₅ , 99%; ¹⁵ N, 99%	2.5
L-Phenylalanine	Phe	¹³ C ₉ , 99%; ¹⁵ N, 99%	2.5
L-Proline	Pro	¹³ C ₅ , 99%; ¹⁵ N, 99%	2.5
L-Serine	Ser	¹³ C ₃ , 99%; ¹⁵ N, 99%	2.5
L-Threonine	Thr	¹³ C ₄ , 97-99%; ¹⁵ N, 97-99%	2.5
L-Tryptophan*	Trp	¹³ C ₁₁ , 99%; ¹⁵ N ₂ , 99%	2.5
L-Tyrosine	Tyr	¹³ C ₉ , 99%; ¹⁵ N, 99%	2.5
L-Valine	Val	¹³ C ₅ , 99%; ¹⁵ N, 99%	2.5

^{*}Compounds absent in MSK-A2.

Companion unlabeled standard mixes and kits may be available; please inquire.




Figure. Application of the MSK-A2 mix to compare the quantitative profiles of cow serum and standard media (see PMID: 28826492 for additional information).

Example References

Gao, J.; Chung, T-S. 2021. Membranes made from nonsolvent-thermally induced phase separation (N-TIPS) for decellularization of blood in dry plasma spot (DPS) applications. Chem Eng Sci, 229. doi.org/10.1016/j.ces.2020.116010.

van Gastel, N.; Spinelli, J.B.; Sharda, A.; et al. 2020. Induction of a timed metabolic collapse to overcome cancer chemoresistance. Cell Metab, 32(3), 391-403.

Røst, L.M.; Thorfinnsdottir, L.B.; Kumar, K.; et al. 2020. Absolute quantification of the central carbon metabolome in eight commonly applied prokaryotic and eukaryotic model systems. Metabolites, 10(2), 74.

Alcock, R.D.; Shaw, G.C.; Tee, N.; et al. 2019. Plasma amino acid concentrations after the ingestion of dairy and collagen proteins, in healthy active males. Front Nutr, 6, 163. Sullivan, M.R; Danai, L.V.; Lewis, C.A.; et al. 2019. Quantification of microenvironmental metabolites in murine cancers reveals determinants of tumor nutrient availability. Elife, 8, e44235.

Chen, W.W.; Freinkman, E.; Sabatini, D.M. 2017. Rapid immunopurification of mitochondria for metabolite profiling and absolute quantification of matrix metabolites. Nat Protoc, 12(10), 2215-2231.

Muir, A.; Danai, L.V.; Gui, D.Y.; et al. 2017. Environmental cystine drives glutamine anaplerosis and sensitizes cancer cells to glutaminase inhibition. Elife, 6, e27713.

MSK-NCAA Mix

Composition

Reconstituting the dried-down NCAA mix with 1 mL solvent (e.g., 50% methanol) results in 2.5 mM concentrations for each compound.

Compound	Abbrev.	Label and Enrichment
β-Alanine	β-Ala	¹³ C ₃ , 98%; ¹⁵ N, 96-99%
L-Azidohomoalanine·HCl	hAHA	1,2,3,4- ¹³ C ₄ , 99%; 2,4- ¹⁵ N ₂ , 98%
L-Citrulline	Cit	1,2,3,4,5- ¹³ C ₅ , 98%
L-Dihydroxyphenylalanine	DOPA	1-13C, ring-13C ₆ , 99%
L-Homoarginine·HCl	Harg	¹³ C ₇ , 98%; ¹⁵ N ₄ , 98%
L-Ornithine·HCl	Orn	¹³ C ₅ , 98%
Sarcosine·HCl	Sar	¹³ C ₃ , 99%; ¹⁵ N, 98%

Usage Specifications

Mix Type	MSK-A2	MSK-CAA	MSK-NCAA
Form	1.2 mL solution	dried	down
Before reconstitution:			
Storage	-5 to 5°C; protect from light		
Recommended retest	2 years from date of manufacture		
Upon reconstitution:			
Storage	N/A	-5 to 5°C; pro	tect from light
Recommended retest	N/A	4 w	reeks

I use several of CIL's metabolomics mixes in my LC-MS analysis. They make quantitative metabolomic work convenient with a single internal standard mix spike and provide the corresponding unlabeled mix for absolute quantitation. The ready mixes save time from the tedious task of making up individual solutions of each analyte for quantitation. I find the different mixes applicable to a variety of analyses in complicated matrices that provide metabolic insight to my studies.

> - Andrew Downey Senior Scientist, Axcella Health Inc. (USA)

Companion unlabeled standard mixes and kits may be available; please inquire.

Chemical purity (CP) is 98% or greater, unless otherwise indicated

NSK-A Mix

Composition

Reconstituting a given NSK-A vial's contents in 1 mL of high-purity solvent (50:50 water:methanol recommended) will produce the tabulated concentrations below. To facilitate complete dissolution, it is recommended to vortex manually for 1 minute then auto-vortex for a minimum of 30 seconds. **Note:** A combined set of NSK-A and NSK-B – the carnitine/acylcarnitine reference standard mix – is also available (NSK-AB).

Compound	Abbrev.	Label and Enrichment	Conc. (µM)
L-Alanine	Ala	2,3,3,3-D ₄ , 98%	500
L-Arginine·HCI	Arg	5-13C, 99%; 4,4,5,5-D ₄ , 95%	500
L-Aspartatic acid	Asp	2,3,3-D ₃ , 98%	500
L-Citruline	Cit	5,5-D ₂ , 98%	500
DL-Glutamatic acid	Glu	2,4,4-D ₃ , 98%	500
Glycine	Gly	2- ¹³ C, 99%; ¹⁵ N, 98%	2500
L-Leucine	Leu	5,5,5-D ₃ , 99%	500
L-Methionine	Met	methyl-D ₃ , 98%	500
L-Ornithine·HCI*	Orn	5,5-D ₂ , 98%	500
L-Phenylalanine	Phe	ring- ¹³ C ₆ , 99%	500
L-Tyrosine	Tyr	ring- ¹³ C ₆ , 99%	500
L-Valine	Val	D ₈ , 98%	500

^{*}NSK-A1 contains Orn 3,3,4,4,5,5,- D_{G} 98% instead of 5,5- D_{2} , 98%. The remaining components and concentrations are equivalent.

Usage Specifications

Criteria	Recommendation
No. of uses	960 samples/vial
Before reconstitution:	
Storage	≤25°C; protect from light
Recommended retest	4 years from date of manufacture
Upon reconstitution:	
Storage	5±3°C in a tightly sealed vial Note: Storing the sealed vial in a second sealed container helps maintain the integrity of the solution.
Recommended retest	4 weeks

We have been using CIL's NSK-A and NSK-B reference standards in our LC-MS/MS method for newbom screening tests since 2001. From that time more than 4,000,000 newborns have been screened. The products maintain good stability after dissolution and the high analysis repeatability enabled us to stabilize their use as internal standards. We appreciate the high-quality products, timely deliveries, and excellent customer relations.

Mariusz Oltarzewski
 Head of Department of Screening Tests and Metabolic Diagnostics
 Institute of Mother and Child (Poland)




Figure. Full scan MS1 spectra of the DBS internal standards, with the derivatized standards shown in A and the nonderivatized in B. The standard mixture consists of stable isotope-labeled amino acids (NSK-A and -T) and carnitine/acylcarnitines (NSK-B and NSK-B-G). Please refer to Thermo Scientific technical note #73398 for method and application details.

Example References

Zhang, S.; Li, X.; Luo, H.; et al. 2020. Role of aromatic amino acids in pathogeneses of diabetic nephropathy in Chinese patients with type 2 diabetes. J Diabetes Complications, 34(10), 107667.

Wang, L.; Liu, D.; Shen, H.; et al. 2020. Analysis of amino acid patterns with nutrition regimens in preterm infants with extrauterine growth retardation. Front Pediatr, 8, 184. Brennenstuhl, H.; Kohlmüller, D.; Gramer, G.; et al. 2020. High throughput newborn screening for aromatic L-amino-acid decarboxylase deficiency by analysis of concentrations of 3-O-methyldopa from dried blood spots. J Inherit Metab Dis, 43(3), 602-610.

Cao, Y-F.; Li, J.; Zhang, Z; et al. 2019. Plasma levels of amino acids related to urea cycle and risk of type 2 diabetes mellitus in Chinese adults. Front Endocrinol, 10, 50. Jing, F.; Hu, X.; Cao, Y.; et al. 2018. Discriminating gastric cancer and gastric ulcer using human plasma amino acid metabolic profile. UBMB Life, 70(6), 553-562.

Xie, X.; Kozak, M. 2020. Simultaneous analysis of amino acids, acylcarnitines, and succinylacetone in dried blood spots for research using nonderivatized and derivatized methods. (Thermo Scientific #73398).

Companion unlabeled standard mixes and kits may be available; please inquire.

Bile Acid Mixes

Bile acids (BAs) are steroid-like compounds that act as a detergent in the breakdown of fats. This family of compounds comprises primary BAs (synthesized in the liver) and secondary BAs (produced in the gut by bacteria). These are essential regulatory compounds that are involved in various metabolic processes (e.g., cholesterol and lipid metabolism) and signaling interactions (e.g., in glucose and energy homeostasis). Investigations into their synthesis/ metabolism, disease linkage, and biomarker potential are examples of the type of research studies being undertaken. To aid further research and development efforts in this space, CIL has formulated stable isotope-labeled (and unlabeled) BA mixes. These dried-down mixes are constructed with the unconjugated BAs in one vial and the conjugated BAs in a second.

Metabolomics Bile Acid Motures for Custing and Superior Custing Superior C

Click on the thumbnail or visit isotope.com/applications/ for

Overview

Catalog No.	Description	No. of Metabolites	Unit Size
MSK-BA1	Bile Acid Standard Mix 1 – Unconjugated	6	1 vial
MSK-BA2	Bile Acid Standard Mix 2 – Conjugated	10	1 vial

Compositions

Reconstituting a given BA mix in 1 mL of solvent (e.g., 50% methanol) yields a concentration of ~100 μM.

Unconjugated BA Mix (MSK-BA1)				
Compound	Abbrev.	Type of BA	Label and Enrichment	~Qty. (µg)
Chenodeoxycholic acid	CDCA	Primary	2,2,4,4-D ₄ , 98%	40
Cholic acid	CA	Primary	2,2,4,4-D ₄ , 98%	41
Deoxycholic acid	DCA	Secondary	2,2,4,4-D ₄ , 98%	40
Lithocholic acid	LCA	Secondary	2,2,4,4-D ₄ , 98%	38
β-Muricholic acid	β-МСА	Primary	2,2,3,4,4-D ₅ , 99%	41
Ursodeoxycholic acid	UDCA	Secondary	2,2,4,4-D ₄ , 98% (CP 95%)	40

Conjugated BA Mix (MSK-BA2)				
Compound	Abbrev.	Type of BA	Label and Enrichment	~Qty. (µg)
Glycochenodeoxycholic acid	GCDCA	Primary	2,2,4,4-D ₄ , 98% (CP 97%)	45
Glycocholic acid	GCA	Primary	2,2,4,4-D ₄ , 98% (CP 96%) contains ~4% water	47
Glycodeoxycholic acid	GDCA	Secondary	2,2,4,4-D ₄ , 98%	45
Glycolithocholic acid	GLCA	Secondary	2,2,4,4-D ₄ , 98%	44
Glycoursodeoxycholic acid	GUDCA	Secondary	2,2,4,4-D ₄ , 98% (CP 97%)	45
Taurochenodeoxycholic acid, sodium salt	TCDCA	Primary	2,2,4,4-D ₄ , 98% (CP 97%)	53
Taurocholic acid, sodium salt	TCA	Primary	2,2,4,4-D ₄ , 98%	54
Taurodeoxycholic acid, sodium salt	TDCA	Secondary	2,2,4,4-D ₄ , 98%	53
Taurolithocholic acid, sodium salt	TLCA	Secondary	2,2,4,4-D ₄ , 98%	51
Tauroursodeoxycholic acid, sodium salt	TUDCA	Secondary	2,2,4,4-D ₄ , 98%	53

Usage Specifications

Criteria	Recommendation
Before reconstitution:	
Storage	-20°C; protect from light
Recommended retest	2 years from date of manufacture
After reconstitution:	
Storage	4°C in a tightly sealed vial Note: Storing the sealed vial in a second sealed container helps maintain the integrity of the solution.
Recommended retest	4 weeks

Companion unlabeled standard mixes and kits may be available; please inquire.




Figure. Chromatographic overview of a combined BA mix (six unconjugates and 10 conjugates) analysis measured by RPLC-MS (Orbitrap ID-X, negative ESI). The m/z displayed are the [M-H]⁻ ions. Procedurally. the labeled/unlabeled BA mixes were reconstituted in 50% methanol before aliquot mixing and MS1 measurement. Note: The labeled and unlabeled BAs coeluted.

The materials that Cambridge Isotope Laboratories (CIL) supply are always of the highest quality in terms of chemical and isotopic purity. Recently, we have embarked on two new major research areas involving the analyses of bile acids (BAs) and per- and polyfluorinated alkylsubstances (PFAS). The standards/mixes from CIL that enable these, and other, projects are precisely what we have come to expect perfection!

- David C. Muddiman, PhD Jacob and Betty Belin Distinguished Professor, Department of Chemistry, North Carolina State University Director of Molecular Education, Technology, and Research Innovation Center (METRIC)

Carnitine/Acylcarnitine Mixes and Sets

Carnitine and acylcarnitines play an essential role in fatty acid metabolism. Metabolism disorders of fatty acid oxidation and several organic acidurias impose major clinical manifestations (e.g., hypoketotic hypoglycemia, skeletal myopathy, liver disease, and/or failure). These are largely attributed to enzymatic deficiencies and can be monitored through carnitine/acylcarnitine measurement. To aid MS/MS screening studies, CIL has formulated stable isotopelabeled (and unlabeled) standard carnitine/acylcarnitine mixes.

Overview

Catalog No.	Description	No. of Metabolites	Unit Size
NSK-B*	Carnitine/Acylcarnitine Standard Mix Set B	8	1 vial, 10 vials
NSK-B-G1	Carnitine/Acylcarnitine Standard Mix Supplement to NSK-B	5	1 vial, 10 vials
NSK-AB	Standard Mix Sets A and B	12 (in A), 8 (in B)	2 × 10 vials



Click on the thumbnail or visit isotope.com/applications/ for

In clinical chemistry-based applications of mass spectrometry, the first lesson the laboratory learns is the requisite nature of stable isotopeenriched standards for quantification of metabolites in biological fluids. In newborn screening of amino acids and acylcarnitines, Cambridge Isotope Laboratories, Inc., set the standard for quantification of these metabolites in dried blood spots. As research and development of the newborn screening analysis by mass spectrometry progressed, it was clear that a half dozen isotope-labeled internal standards would not be adequate for the analysis of an amino acid and acylcarnitine profile, together comprising a range of 500 separate mass units and more than 30 important metabolites, most of which require accurate quantification. When screening began to expand beyond research, it was clear that weighing out small quantities of individual standards would reduce accuracy and introduce unnecessary error. Therefore, together, we set out to develop sets of standards for amino acids and acylcarnitine analysis that would enable quantification. We started this development more than 20 years ago adding, changing and improving these standards. CIL, together with the early developers of tandem mass spectrometrybased newborn screening, set the standard by which all other laboratories follow. CIL's commitment to supporting the metabolic and newborn screening community is exceptional. It is our good fortune in the clinical chemistry and mass spectrometry community to have CIL as part of our laboratory solutions.

> - Donald H. Chace, PhD MSFS FACB Medolac Laboratories (USA)

Companion unlabeled standard mixes and kits may be available; please inquire.

Chemical purity (CP) is 98% or greater, unless otherwise indicated

Compositions

Reconstituting a given NSK-B or NSK-B-G1 vial's contents in 1 mL of highly pure methanol will produce the tabulated concentrations below. To facilitate complete dissolution, it is recommended to vortex manually for 1 minute then auto-vortex for a minimum of 30 seconds. **Note:** A combined set of NSK-B and NSK-A – the amino acid reference standard mix – is also available (NSK-AB).

NSK-B				
Compound	Abbrev.	Label and Enrichment	Conc. (µM)	
L-Carnitine	CO	trimethyl-D ₉ , 98%	152	
O-Acetyl-L-carnitine·HCl	C2	N-methyl-D ₃ , 98%	38	
O-Propionyl-L-carnitine·HCl	C3	N-methyl-D ₃ , 98%	7.6	
O-Butyryl-L-carnitine·HCl	C4	N-methyl-D ₃ , 98%	7.6	
O-Isovaleryl-L-carnitine·HCl	C5	N,N,N-trimethyl-D ₉ , 98%	7.6	
O-Octanoyl-L-carnitine·HCl	C8	N-methyl-D ₃ , 98%	7.6	
O-Myristoyl-L-carnitine·HCl	C14	N,N,N-trimethyl-D ₉ , 98%	7.6	
O-Palmitoyl-L-carnitine·HCl	C16	N-methyl-D ₃ , 98%	15.2	

NSK-B-G1					
Compound Abbi		Label and Enrichment	Conc. (µM)		
O-Glutaryl-L-carnitine⋅ClO ₄	C5-DC	<i>N</i> -methyl-D ₃ , 98% (CP 97%)	15.2		
3-Hydroxyisovaleryl-L-carnitine·ClO ₄	C5-OH	N-methyl-D ₃ , 98%	7.6		
O-Dodecanoyl-L-carnitine·HCl	C12	N,N,N-trimethyl-D ₉ , 98%	7.6		
O-3-DL-Hydroxypalmitoyl-L-carnitine·ClO ₄	C16-OH	N-methyl-D ₃ , 98%	15.2		
O-Octadecanoyl-L-carnitine·HCl	C18	N-methyl-D₃, 98%	15.2		

Usage Specifications

Mix Type	NSK-B	NSK-B-G1	
No. of uses	960 samples/vial		
Before reconstitution:			
Storage	≤8°C; protec	ct from light	
Recommended retest	1 year from date of manufacture 2 years from date of manufacture		
After reconstitution:			
Storage	5±3°C in a tightly sealed vial	5±3°C in a tightly sealed vial	
	Note: Storing the sealed vial in a second sealed container helps maintain the integrity of the solution		
Recommended retest	4 weeks		




Figure. Full scan MS1 spectra of the DBS internal standards, with the derivatized standards shown in **A** and the nonderivatized in **B**. The standard mixture consists of stable isotope-labeled amino acids (NSK-A and -T) and carnitine/acylcarnitines (NSK-B and NSK-B-G). Please refer to Thermo Scientific technical note #73398 for method and application details.

Example References

Brailova, M.; Clerfond, G.; Trésorier, R.; et al. **2020**. Inherited metabolic diseases and cardiac pathology in adults: diagnosis and prevalence in a cardiometabo study. *J Clin Med*, *9*(3), 694.

Varkey, A.; Devi, S.; Mukhopadhyay, A.; et al. **2020**. Metabolome and microbiome alterations related to short-term feeding of a micronutrient-fortified, high-quality legume protein-based food product to stunted school age children: a randomized controlled pilot trial. *Clin Nutr, 39(11), 3251-3261*.

Brennenstuhl, H.; Kohlmüller, D.; Gramer, G.; et al. **2020**. High throughput newborn screening for aromatic t-amino-acid decarboxylase deficiency by analysis of concentrations of 3-O-methyldopa from dried blood spots. *J Inherit Metab Dis*, *43*(3), 602-610.

Cao, B.; Wang, D.; Pan, Z.; et al. **2019**. Characterizing acylcarnitine biosignatures for schizophrenia: a longitudinal pre- and post-treatment study. *Transl Psychiatry*, *9*(1), 19. Puskarich, M.A.; Evans, C.R.; Karnovsky, A.; et al. **2018**. Septic shock nonsurvivors have persistently elevated acylcarnitines following carnitine supplementation. *Shock*, *49*(*4*), 412-419.

Technical Note

Xie, X.; Kozak, M. 2020. Simultaneous analysis of amino acids, acylcarnitines, and succinylacetone in dried blood spots for research using nonderivatized and derivatized methods. (Thermo Scientific #73398).

Companion unlabeled standard mixes and kits may be available; please inquire.

Chemical purity (CP) is 98% or greater, unless otherwise indicated.

Fatty Acid/Lipid Mixes

Fatty acids and lipids are important biological compounds that are essential to the regulation and control of cellular functions and metabolic pathways. These biomolecules are also tied to the energetic balance of an organism. Their qualitative/quantitative analysis has emerged to better understand the underlying pathophysiology, as well as to identify new biomarkers or diagnose existing ones. To aid such initiatives, CIL offers an array of mixed fatty acids and triglycerides. These dried-down mixes are uniformly labeled and available in different forms (i.e., free acid, methyl ester) as research-grade material.

Click on the thumbnail or visit isotope.com/applications/ for

Overview

Catalog No.	Description	No. of Metabolites	Unit Size
CLM-8455	Fatty Acid Mix (U-13C, 98%)	5-10	0.25 g, 1 g
DLM-8572	Fatty Acid Mix (U-D, 96-98%)	5-10	Please inquire
CDLM-8376	Fatty Acid Mix (U-13C, 98%; U-D, 97%)	5-10	0.25 g, 1 g
CLM-8381	Fatty Acid, Methyl Ester Mix (U-13C, 98%) (unlabeled terminal ester) CP 95%	5-10	0.25 g, 1 g
DLM-2497	Fatty Acid, Methyl Ester Mix (U-D, 96-98%)	5-10	Please inquire
DLM-8375	Triglyceride Mix (U-D, 97%)	5-10	0.25 g, 1 g

Compositions

The mixed fatty acid and triglyceride products are derived from the Agmenelum quadriplicatum algal source. Please inquire for details as the mix compositions can vary by lot.

Usage Specifications

Min Tuno	Catalag Na	Before Reconstitution		
Mix Type	Catalog No.	Storage	Recommended Retest	
Fatty acids	CLM-8455 DLM-8572 CDLM-8376	-5 to 5°C; protect from light	2 years	
Fatty acid, methyl esters	CLM-8381 DLM-2497	-5 to 5°C; protect from light	5 years	
Triglycerides	DLM-8375	-5 to 5°C; protect from light	5 years	

Example References

Hernandez-Saavedra, D.; Sanders, L.; Freeman, S.; et al. 2020. Stable isotope metabolomics of pulmonary artery smooth muscle and endothelial cells in pulmonary hypertension and with TGF-beta treatment. Sci Rep., 10(1), 413.

He, C.; Weston, T.A.; Junq, R.S.; et al. 2018. NanoSIMS analysis of intravascular lipolysis and lipid movement across capillaries and into cardiomyocytes. Cell Metab, 27(5), 1055-1066.

Nucleic Acid Sets

Nucleic acids are necessary biomolecules of living systems, being fundamentally important to a multitude of cellular processes. Its basic building blocks are nucleobases (e.g., adenine, cytosine), nucleosides (e.g., adenosine, guanosine), and nucleotides (e.g., ATP, CDP). The qualification/quantification of these compounds is conducted for a number of purposes. This includes the screening of metabolic errors and evaluating the efficacy of drug treatments (be it anticancer, antiviral, or immunosuppressive), among other target areas. CIL offers a variety of stable isotope-labeled nucleic acid sets to help advance MS- and/or NMR-based research. The sets consist of deoxyribonucleoside monophosphates (dNMPs), deoxyribonucleoside triphosphates (dNTPs), and ribonucleoside triphosphate (rNTPs) in different unit configurations.

Overview

Catalog No.	Description	No. of Metabolites	Unit Size
CNLM-7871-SL	Set of 4 2'-Deoxyribonucleoside 5'-Monophosphates (U-13C, 98%; U-15N, 98%) (Li salts, in solution) CP 90%	4	10 mg
NLM-7512-SL	Set of 4 2'-Deoxyribonucleoside 5'-Triphosphates (U-15N, 98%) (Li salts, in solution) CP 90%	4	10 mg, 50 mg
DLM-7511-SL	Set of 4 2'-Deoxyribonucleoside 5'-Triphosphates (U-D, 98%) (Li salts, in solution) CP 90%	4	10 mg, 50 mg
CNLM-7513-SL	Set of 4 2'-Deoxyribonucleoside 5'-Triphosphates (U-13C, 98%; U-15N, 98%) (Li salts, in solution) CP 90%	4	10 mg, 50 mg
NLM-7519-SL	Set of 4 Ribonucleoside 5'-Triphosphates (U-15N, 98%) (Li salts, in solution) CP 90%	4	10 mg, 50 mg
DLM-7518-SL	Set of 4 Ribonucleoside 5'-Triphosphates (U-D, 98%) (Li salts, in solution) CP 90%	4	10 mg, 50 mg
CNLM-7503-SL	Set of 4 Ribonucleoside 5'-Triphosphates (U-13C, 98%; U-15N, 98%) (Li salts, in solution) CP 90%	4	10 mg, 50 mg
NLM-7519-CA	Set of 4 Ribonucleoside 5'-Triphosphates (U-15N; 98%) (NH ₄ salts, in solution) CP 90%	4	4 × 100 μmol
DLM-7518-CA	Set of 4 Ribonucleoside 5'-Triphosphates (U-D, 98%) (NH ₄ salts, in solution) CP 90%	4	4 × 20 μmol 4 × 50 μmol 4 × 100 μmol
CNLM-7503-CA	Set of 4 Ribonucleoside 5'-Triphosphates (U-13C, U-15N; 98-99%) (NH ₄ salts, in solution) CP 90%	4	4 × 20 μmol 4 × 50 μmol 4 × 100 μmol

dNMP

Composition and Specification

Catalog No.	dNMP	Conc. (mM)
CNLM-7871-SL	AMP, CMP, GMP, TMP	80

- Supplied on dry ice as a solution of 5 mM Tris-HCl (pH 7.5) in water.
- Store below -20°C; protect from light.

dNTPs

Compositions and Specifications

Catalog No.	dNTP	Conc. (mM)
NLM-7512-SL		
DLM-7511-SL	datp, dctp, dgtp, dttp	50-100
CNLM-7513-SL		

- Supplied on dry ice as a solution of 5 mM Tris-HCl (pH 7.5) in water (for NLM-7512-SL and CNLM-7513-SL) or D₂O (for DLM-7511-SL)
- Store below -20°C; protect from light.

Example References

Liu, B.; Winkler, F.; Herde, M.; et al. 2019. A link between deoxyribonucleotide metabolites and embryonic cell-cycle control. Curr Biol, 29(7), 1187-1192. Song, Y.; Marmion, R.A.; Park, J.O.; et al. 2017. Dynamic control of dNTP synthesis in early embryos. Dev Cell, 42(3), 301-308.

rNTPs

Compositions and Specifications

Catalog No.	rNTP	Conc. (mM)	
NLM-7519-SL			
DLM-7518-SL	ratp, rctp, rgtp, rttp	50-100	
CNLM-7503-SL			

- Supplied on dry ice as a solution of 5 mM Tris-HCl (pH 7.5) in water (for NLM-7519-SL and CNLM-7503-SL) or D₂O (for DLM-7518-SL)
- Store below -20°C; protect from light.

Catalog No.	rNTP	Conc. (mM)
NLM-7519-CA		
DLM-7518-CA	ratp, rctp, rgtp, rutp	20, 50, and/or 100
CNLM-7503-CA		

- Supplied on dry ice as a solution (e.g., 100 mM contains 100 µmol of each rNTP in 1 mL water, pH ~7.5).
- Store at -20°C; protect from light. Recommended retest is 5 years from date of manufacture.

Organic Acid Mixes

Organic acids (OAs) are byproducts of AAs and are intermediates in various biochemical pathways (e.g., glycolysis, TCA cycle). These compounds have emerged as important regulators of host immune response and transcriptional regulation. To help facilitate this analysis in various research applications, CIL has formulated a collection of OA mixes in their stable isotope-labeled and unlabeled forms. These dried-down mixes are available as a master mix (comprises 33 OAs), as well as in subsets (with 33 OAs spread across seven subclasses).

Overview

Catalog No.	Description	Mix Type	No. of Metabolites	Unit Size
MSK-OA	Organic Acid Mix	Master	33	1 vial
MSK-OA1	Keto Acid Mix	Subset 1	5	1 vial
MSK-OA2	Diacid Mix	Subset 2	8	1 vial
MSK-OA3	Monoacid Mix	Subset 3	3	1 vial
MSK-OA4	Hydroxy Acid Mix	Subset 4	5	1 vial
MSK-OA5	Aromatic Acid Mix	Subset 5	5	1 vial
MSK-OA6	Other Acid Mix	Subset 6	2	1 vial
MSK-OA7	Other Organic Acid Mix	Subset 7	5	1 vial



isotope.com/applications/ for

Compositions

Composition of stable isotope-labeled organic acid mixtures. Reconstituting these dried-down mixes (either complete or subsets) in 1 mL of solvent (e.g., water) yields an equimolar concentration of 250 µM.

MSK-OA-1 Labeled Organic Acid Mix (single vial of all 33 compounds)

Subsets

MSK-OA1-1 Labeled Keto Acid Mix	MSK-OA4-1 Labeled Hydroxy Acid Mix
α -Ketobutyric acid, sodium salt ($^{13}C_4$, 98%)	Glycolic acid, sodium salt (1,2-13C ₂ , 99%)
α -Ketoglutaric acid, disodium salt (1,2,3,4- 13 C ₄ , 99%) CP 97%	Malic acid, disodium salt· H_2O ($^{13}C_4$, 99%)
α -Ketoisocaproic acid, sodium salt ($^{13}C_6$, 99%)	Sodium D-3-hydroxybutyrate (13C ₄ , 99%) CP 97%
α -Ketoisovaleric acid, sodium salt ($^{13}C_5$, 98%)	Sodium L-lactate (13C ₃ , 98%) 20% w/w in H ₂ O
Sodium pyruvate (13C ₃ , 99%)	DL 2-Hydroxyglutarate, disodium salt (13C ₅ , 99%)
MSK-OA2-1 Labeled Diacid Mix	MSK-OA5-1 Labeled Aromatic Acid Mix
Adipic acid, disodium salt (13C ₆ , 99%)	Hippuric acid (ring- ¹³ C ₆ , 99%)
Fumaric acid, disodium salt (13C ₄ , 99%)	Homovanillic acid, sodium salt (1,2-13C ₂ , 98-99%)
Maleic acid, disodium salt·H ₂ O (¹³C ₄ , 99%)	Phthalic acid, disodium salt (13C ₄ , 99%)
Malonic acid, disodium salt (13C ₃ , 99%)	Sodium benzoate (ring-13C ₆ , 99%)
Methylmalonic acid, disodium salt (13C ₄ , 99%)	DL-Vanilmandelic acid (ring- ¹³ C ₆ , 99%)
Oxalic acid, disodium salt (1,2-13C ₂ , 99%)	MSK-OA6-1 Labeled Other Acid Mix
Sodium isobutyrate (13C ₄ , 99%)	L-Ascorbic acid (¹³ C ₆ , 99%)
Succinic acid, disodium salt (13C ₄ , 99%)	Sodium D-gluconate (13C ₆ , 99%)
MSK-OA3-1 Labeled Monoacid Mix	MSK-OA7-1 Labeled Other Organic Acid Mix
Sodium acetate (1,2-13C ₂ , 99%)	trans-Aconitic acid (2,4,4'-13C ₃ , 99%) CP 95%
Sodium butyrate (13C ₄ , 99%)	Creatine (15N ₃ , 98%)
Sodium propionate (13C ₃ , 99%)	Orotic acid, sodium salt (15N ₂ , 98%)
	Trisodium citrate, hemihydrate (1,5,6-carboxy-13C ₃ , 99%)
	Uric acid, sodium salt (15N ₂ , 98%) CP 95%

Usage Specifications

Criteria	Recommendation	
Before reconstitution:		
Storage	ambient temperature; protect from light and moisture	
Recommended retest	2 years from date of manufacture	

Example Results




Figure. XICs of example metabolites measured in the labeled OA master mix by HILIC-MS (Orbitrap ID-X, negative ESI). Shown are a collection of OAs from a variety of subset classes. The m/z displayed are the [M-H] ions. Procedurally, the dried-down master mix was reconstituted in 1 mL of 50% methanol then an aliquot 10-fold diluted before MS1 analysis.

Application Note

Petucci, C.*; Percy, A.J.*; Zelenin, A.; Gardell, S.J.; Backiel, K. 2017. Organic acid quantitation in mouse muscle by ion chromatography-mass spectrometry with isotopically labeled standards. (Cambridge Isotope Laboratories, Inc. Application Note 47)

^{*}These authors contributed equally to this application note.

Other Mixes

A number of additional stable isotope-labeled mixes are also available at CIL (see overview). Included in this collection is a lysophosphatidylcholine (LPC or lysoPC) mix (NSK-LPC) and MS/MS screening mixes for basic and translational MS research. The details of these mix offerings are outlined below.

Overview

Catalog No.	Description	No. of Metabolites	Unit Size
NSK-LPC	Lysophosphatidylcholine Mix	4	1 vial
NSK-NI	Acid Sphingomyelinase Substrate and Internal Standard Mix	1 (S + IS)	1 vial
NSK-KR	Galactocerebrosidase Substrate and Internal Standard Mix	1 (S + IS)	1 vial
NSK-FA	α-Galactosidase Substrate and Internal Standard Mix	1 (S + IS)	1 vial
NSK-GA	Glucocerebrosidase Substrate and Internal Standard Mix	1 (S + IS)	1 vial
NSK-MP	α-L-Iduronidase Substrate and Internal Standard Mix	1 (S + IS)	1 vial
NSK-PO	Acid α-Glucosidase Substrate and Internal Standard Mix	1 (S + IS)	1 vial



Click on the thumbnail or visit isotope.com/applications/ for

NSK-LPC Mix

Composition

Reconstituting a given LysoPC mix (labeled, NSK-LPC; unlabeled, NSK-LPC-US) in 1 mL of high-purity solvent (e.g., 95:5 v/v methanol:water) will produce the concentrations specified below. To facilitate complete dissolution, it is recommended to sonicate the vial for 3 minutes then auto-vortex for a minimum of 10 seconds.

Standard (Abbreviation)	Label and Enrichment	Conc. (µM)	Structure
Lysophosphatidylcholine 20:0 (LysoPC 20:0)	eicosanoyl-12,12,13,13-D ₄ , 98%	5.5	CH ₃ (CH ₂) ₆ (CD ₂) ₂ (CH ₂) ₁₀ O O O O O O O O O O O O O O O O O O O
Lysophosphatidylcholine 22:0 (LysoPC 22:0)	docosanoyl-1,2,3,4,5,6- ¹³ C ₆ , 99%	5.5	CH ₃ (CH ₂) ₁₅ (CH ₂) ₅ OH ₂ OH ₃ CH ₃
Lysophosphatidylcholine 24:0 (LysoPC 24:0)	tetracosanoyl-1,2,3,4,5,6- ¹³ C ₆ , 99%	5.5	CH ₃ (CH ₂) ₁₇ (CH ₂) ₅ OH ₂ OH ₃ CH ₃ CH ₃
Lysophosphatidylcholine 26:0 (LysoPC 26:0)	hexacosanoyl-1,2,3,4,5,6- ¹³ C ₆ , 99%	5.5	CH ₃ (CH ₂) ₁₉ (CH ₂) ₅ OH OH CH ₃ CH ₃

Usage Specifications

Criteria	Recommendation		
No. of uses	~765 samples/vial		
Before reconstitution:		After reconstitution:	
Storage	≤20°C; protect from light	Storage	5±3°C or -20±5°C
Recommended retest	1 year from date of manufacture	Recommended retest	6 weeks

Testing 1,200 to 2,000 newborn specimens every day is a challenging job. In order to ensure high data quality and proper efficacy of an MS-based screening assay (1st and 2nd tier), we utilize CIL's lyophilized LysoPC mix. This 4-plex cocktail helps reduce significantly our laboratory day-to-day variation and preparation time. We also use the LysoPC mixes (both unlabeled and isotope-labeled) in the tuning of our mass spectrometers and in determining response factors. This ensures that the concentrations rendered for each LysoPC is accurately obtained and can be cross-checked from platform-10-platform.

 Kuldeep Dhillon, Research Scientist Supervisor I California Department of Public Health (USA)x

Companion unlabeled standard mixes and kits may be available; please inquire.

Chemical purity (CP) is 98% or greater, unless otherwise indicated.




Figure. XICs of four labeled LysoPCs from the NSK-LPC mix measured by LC-MS/MS (positive ESI, Xevo-TQS). These serve as individual calibrators for their unlabeled counterparts in sample screening experiments. Data courtesy of CA Department of Public Health.

MS/MS Screening Mixes

α-Galactosidase Substrate and Internal Standard (NSK-FA-1)

Each vial contains the following compounds at a molar ratio of 500:1.

Substrate		Internal Standard	
(6-Benzoylamino-hexyl)-{2-[4-(3,4,5-trihydroxy-6-hydroxymethyl-tetrahydro-pyran-2-yloxy)-phenylcarbamoyl]-ethyl}-carbamic acid <i>tert</i> -butyl ester		(6-D ₅ -Benzoylamino-hexy acid <i>tert</i> -butyl ester	yl)-[2-(4-hydroxy-phenyl-carbamoyl)-ethyl]-carbamic
C ₃₃ H ₄₇ N ₃ O ₁₀	MW: 645.7 Da	C ₂₇ H ₃₂ N ₃ O ₅ D ₅	MW: 488.5 Da
HO OH HO CH3 OCH3 OCH3		но	$- \bigcup_{\mathbf{H}} \bigvee_{\mathbf{O}} \bigvee_{\mathbf{O}} \bigcup_{\mathbf{C}} \bigcup_{\mathbf{C}} \bigcup_{\mathbf{H}_3}^{\mathbf{D}} \bigvee_{\mathbf{O}} \bigcup_{\mathbf{D}} \bigcup_{\mathbf{C}} \bigcup_{\mathbf{D}} \bigcup_{\mathbf{C}} \bigcup_{\mathbf{C}}$

Glucocerebrosidase Substrate and Internal Standard (NSK-GA-1)

Each vial contains the following compounds at a molar ratio of 50:1.

Substrate		Internal Standard		
D-Glucosyl-β1-1'- <i>N</i> -dodecanoyl-D- <i>erythro</i> -sphingosine [C12-glucocerebroside]		N-Myristoyl-D-erythro-s	N-Myristoyl-D- <i>erythro</i> -sphingosine [C14-ceramide]	
C ₃₆ H ₆₉ NO ₈	MW: 643.9 Da	C ₃₂ H ₆₃ NO ₃	MW: 509.8 Da	
HO OH OH OH OCH3			HO HN CH ₃	

Galactocerebrosidase Substrate and Internal Standard (NSK-KR-1)

Each vial contains the following compounds at a molar ratio of 150:1.

Substrate		Internal Standa	ard
D-Galactosyl-β1-1'-octanoyl-D- <i>erythro</i> -sphingosine [C8-galactosylceramide]		N-Decanoyl-D-e	rythro-sphingosine [C10-ceramide]
C ₃₂ H ₆₁ NO ₈	MW: 587.8 Da	C ₂₈ H ₅₅ NO ₃	MW: 453.7 Da
но	OH OHN CH3		$HO \longrightarrow CH_3$

α-L-Iduronidase Substrate and Internal Standard (NSK-MP-1)

Each vial contains the following compounds at a molar ratio of 150:1.

Substrate		Internal Standard	
(7-(1-Iduronic acid)-oxycoumarin-4-methylamine-(5'-N-boc-aminopentanoyl)-amide)		(7-Hydroxycoumarin-4-methylamine- (4'-N-boc-aminobutanoyl)-amide)
C ₂₆ H ₃₄ N ₂ O ₁₂	MW: 566.55 Da	$C_{19}H_{24}N_2O_6$	MW: 376.4 Da
	HO ₂ C C C C C C C C C C C C C C C C C C C	O HN HO	CH ₅ CH ₅ CH ₅

Companion unlabeled standard mixes and kits may be available; please inquire.

Chemical purity (CP) is 98% or greater, unless otherwise indicated

Acid Sphingomyelinase Substrate and Internal Standard (NSK-NI-1)

Each vial contains the following compounds at a molar ratio of 50:1.

Substrate		Internal Standard	
N-Hexanoyl-D-erythro-sphingosylphosphorylcholine [C6-sphingomyelin]		N-Butyroyl-D-erythro	p-sphingosine [C4-ceramide]
C ₂₉ H ₅₉ N ₂ O ₆ P	MW: 562.8 Da	Oa C ₂₂ H ₄₃ NO ₃ MW: 369.6 Da	
CH H ₃ C-N [*] H ₃ C	OH OCH3		OH HN CH ₃

Acid α-Glucosidase Substrate and Internal Standard (NSK-PO-1)

Each vial contains the following compounds at a molar ratio of 100:1.

Substrate	Internal Standard	
(7-Benzoylamino-heptyl)-{2-[4-(3,4,5-trihydroxy-6-hydroxymethyl-tetrahydro-	(7-D ₅ -Benzoylamino-heptyl)-[2-(4-hydroxy-phenyl-carbamoyl)-ethyl]-carbamic	
pyran-2-yloxy)-phenylcarbamoyl]-ethyl}-carbamic acid tert-butyl ester	acid tert-butyl ester	
C ₃₄ H ₄₉ N ₃ O ₁₀ MW: 659.8 Da	C ₂₈ H ₃₄ N ₃ O ₅ D ₅ MW: 502.7 Da	
$\begin{array}{c} OH \\ HO \\ OH \\ OH \\ OH \\ OH \\ OH \\ OH $	$HO- \underbrace{ \begin{array}{c} \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\ \\$	

Usage Specifications

Criteria	Recommendation	
Use	~600 samples/vial	
Before reconstitution:		
Storage	-20°; protect from light	
Recommended retest	2 years from date of manufacture	
After reconstitution:		
Storage	5±3°C or -20±5°C	
Recommended retest	4 weeks	

Example References

Ribas, G.; De Mari, J.F.; Civallero, G.; et al. 2017. Validation of a multiplex tandem mass spectrometry method for the detection of selected lysosomal storage diseases in dried blood spots. JIMES, 5, 1-7.

Tortorelli, S.; Turgeon, C.T.; Gavrilov, D.K.; et al. 2016. Simultaneous testing for 6 lysosomal storage disorders and x-adrenoleukodystrophy in dried blood spots by tandem mass spectrometry. Clin Chem, 62(9), 1248-1254.

Cho, S.E.; Kwak, J.R.; Lee, H.; et al. 2016. Triplex tandem mass spectrometry assays for the screening of 3 lysosomal storage disorders in a Korean population. Clin Chim Acta, 454, 20-27.

Steroid Mixes and Sets

Steroids play vital roles in the regulation of a diverse array of cellular functions and physiological processes. These pertain to development, reproduction, homeostasis, and metabolism, among others. Accurate quantification of this compound class is essential for basic and translation research. To aid MS-based research endeavors in this space, CIL is pleased to offer a few different types of stable isotope-labeled steroid mixes.

Overview

Catalog No.	Description	No. of Metabolites	Unit Size
NSK-S	Steroid Mix Set S	5	1 vial, 10 vials
NSK-S-40X	Steroid Mix Set S (40X)	5	1 vial, 10 vials
NSK-S-EXP	Expanded Steroid Mix Set S	9	1 vial, 10 vials



Click on the thumbnail or visit isotope.com/applications/ for

NSK-S Mix

Composition

Reconstituting a given vial's contents in 1 mL of highly pure methanol will produce the concentrations listed in the table below. To facilitate complete dissolution, it is recommended to vortex manually for 1 minute then auto-vortex for a minimum of 30 minutes. Note: NSK-S-40X has 40x the specified concentrations.

Compound	Abbrev.	Label and Enrichment	Conc. (nM)
4-Androstene-3,17-dione	A4	2,2,4,6,6,16,16-D ₇ , 97%	20
Cortisol	F	9,11,12,12-D ₄ , 98%	100
11-Deoxycortisol	11-S	2,2,4,6,6-D ₅ , 98%	20
21-Deoxycortisol	21-S	2,2,4,6,6,21,21,21-D ₈ , 97%	20
17α-Hydroxyprogesterone	17-OHP	2,2,4,6,6,21,21,21-D ₈ , 98%	20

Usage Specifications

Criteria	Recommendation	
No. of uses	48 samples/vial	
Before reconstitution:		
Storage	-5 to 5°C; protect from light	
Recommended retest	5 years from date of manufacture	
Upon reconstitution:		
Storage	Store in a tightly sealed vial at 5±3°C. To maintain the integrity of the solution, store the sealed vials in a second sealed container.	
Recommended retest	4 weeks	

Example References

Gervasoni, J.; Schiattarella, A.; Primiano, A.; et al. 2016. Simultaneous quantification of 17-hydroxyprogesterone, androstenedione, testosterone and cortisol in human serum by LC-MS/MS using TurboFlow online sample extraction. Clin Biochem, 49(13-14), 998-1003.

Hicks, R.A.; Yee, J.K.; Mao, C.S.; et al. 2014. Precursor-to-product ratios reflect biochemical phenotype in congenital adrenal hyperplasia. Metabolomics, 10(1), 123-131.

NSK-S-EXP Mix

Composition

Reconstituting a given vial's contents in 1 mL of highly pure methanol will produce the concentrations listed in the table below. To facilitate complete dissolution, it is recommended to vortex manually for 1 minute then auto-vortex for a minimum of 30 minutes.

Compound	Abbrev.	Label and Enrichment	MW (Da)	Conc. (µM)
Aldosterone	А	D ₇ , 98%	367.49	0.52
4-Androstene-3,17-dione	A4	2,2,4,6,6-D ₅ , 98%	291.44	0.12
Corticosterone	В	2,2,4,6,6,17α,21,21-D ₈ , 97-98%	354.51	1.58
Cortisol	F	9,12,12-D ₃ , 98%	365.48	2.57
Dehydroepiandrosterone sulfate- sodium salt-2H ₂ O	DHEAS	2,2,3,4,4,6-D ₆ , 95%	432.54	21.69
11-Deoxycortisol	11-S	2,2,4,6,6-D ₅ , 98% (CP 97%)	351.49	0.54
17-α-Hydroxyprogesterone	17-OHP	2,2,4,6,6,21,21,21-D ₈ , 98%	338.51	0.27
Progesterone	Р	2,2,4,6,6,17α,21,21,21-D ₉ , 98%	323.52	0.14
Testosterone	Т	2,2,4,6,6-D ₅ , 98%	293.46	0.12

Usage Specifications

Criteria	Recommendation	
Before reconstitution:		
Storage	-5 to 5°C; protect from light	
Recommended retest	1 year from date of manufacture	

Kits

To help facilitate product application, a number of the CIL mix offerings are accompanied with a user manual. The kit manual provides a guide to the preparation and processing possibilities for product implementation, as well as supplying further resources for additional user reference. The majority of the kits supply the supporting documents through a QR code (via product label on kit materials box); the exceptions are the PeptiQuant™ Plus kits (documents on USB) and the IROA kits (documents delivered in email). The accompanying documents minimally include a user manual, certificate of analysis (CoA), safety data sheet (SDS), and product flyer. Overall, the mixes are amenable for use in quality control and qualification/quantification exercises in the MS 'omics space using targeted, semitargeted, or untargeted LC-MS methodologies. Note that the individual mixes may be obtained separately without the manual. Please inquire.

Metabolomic QC Kits

To ensure high-quality metabolomics results, the method and instrument platform must be qualified as being fit-for-purpose. This involves testing for losses or errors in the analytical workflow. To aid such performance assessments in MS metabolomics and potentially other applications (e.g., qualification, quantification), CIL offers two types of QC kits: MSK-QC-KIT and MSK-QReSS-KIT (see details below). These kits comprise two vials of dried-down metabolite mixes (see composition tables below) and have been qualified on both low- and high-resolution mass spectrometers operated under an untargeted/targeted metabolomics regimen (see example results below).

Overview

Catalog No.	Description	No. of Metabolites	Kit Contents
MSK-QC-KIT	Metabolomics QC Kit	5 (in vial 1) 9 (in vial 2)	2 vials of ¹³ C-labeled metabolites document package (user manual, CoA, SDS, product flyer)
MSK-QReSS-KIT*	Metabolomics QReSS Kit	12 (in vial 1) 6 (in vial 2)	 2 vials of isotope labeled metabolites document package (user manual, CoA, SDS, product flyer)

^{*}QReSS™ stands for Quantification, <u>Re</u>tention, and <u>S</u>ystem <u>S</u>uitability. This kit was developed and matrix tested in collaboration with SCIEX.

MSK-QC Kit

Mix Compositions

Reconstituting a given QC mix vial in 1 mL of solvent (e.g., 0.1% FA in 5% ACN) will yield the specified compound concentrations.

Compound	Abbrev.	Label and Enrichment	Conc. (µg/mL)	Vial No.
L-Alanine	Ala	¹³ C ₃ , 99%	4	1
L-Leucine	Leu	¹³ C ₆ , 99%	4	1
L-Phenylalanine	Phe	¹³ C ₆ , 99%	4	1
L-Tryptophan	Trp	¹³ C ₁₁ , 99%	40	1
L-Tyrosine	Tyr	¹³ C ₆ , 99%	4	1
D-Glucose	Glc	¹³ C ₆ , 99%	4	2
D-Sucrose	Suc	¹³ C ₆ , 98%	4	2
Caffeine	CAF	¹³ C ₃ , 99%	4	2
Stearic acid, sodium salt	18:0	¹³ C ₁₈ , 98%	0.4	2
Sodium octanoate	8:0	¹³ C ₈ , 99%	4	2
Sodium propionate	PA	¹³ C ₃ , 99%	4	2
Sodium benzoate	BZA	¹³ C ₆ , 99%	4	2
Sodium citrate	CA	¹³ C ₃ , 99%	4	2
Succinic acid, disodium salt	SA	¹³ C ₄ , 99%	4	2

Usage Specifications

Criteria	Recommendation	
Before reconstitution:		
Storage	ambient temperature	
Recommended retest	2 years from date of manufacture	





Click on the thumbnails or visit isotope.com/applications/ for more information

Companion unlabeled standard mixes and kits may be available; please inquire.




Figure. Representative XICs of a subset of kit metabolites – vial 1 in A and vial 2 in B – measured in human plasma by RPLC-MS (negative ESI, Q Exactive). Note that all isotopically labeled metabolites in the mixes co-eluted with their endogenous analytes in the plasma analyses and their m/z were predominantly [M-H].

I have used products from Cambridge Isotope Laboratories (CIL) for well over 15 years because of the high quality provided. I have collaborated with CIL to develop a new product that would enable improved quality control in MS metabolomics. This process was exciting and engaging. Through collaboration, I found CIL to be a cutting-edge vendor that seeks opportunities to enable scientific discovery and data quality. They seek advice from clients to learn how to better serve them. I always look forward to talking with CIL at conferences and often seek out their booth to visit, not just to talk science, but to also see how their team is doing.

– Timothy J. Garrett, PhD

Associate Professor in Department of Pathology,
Immunology, and Laboratory Medicine, University of Florida (USA)

MSK-QReSS Kit

Compositions

Reconstituting a given QReSS mix vial in 1 mL of solvent (e.g., 50% methanol) will yield the specified compound concentrations.

Compound	Abbrev. or Alt. Name	Label and Enrichment	Conc. (µg/mL)	Vial No.
L-Alanine	Ala	¹³ C ₃ , 99%; ¹⁵ N, 99%	100	1
1,4-Butanediamine·2HCl	putrescine	¹³ C ₄ , 99%	10	1
Creatinine	Crn	N-methyl-D ₃ , 98%	100	1
Ethanolamine·HCl	ETA	1,1,2,2-D ₄ , 98%	10	1
Guanosine·2H ₂ O	Guo	¹⁵ N ₅ , 96-98%	2	1
Hypoxanthine	HPX	¹³ C ₅ , 99%	10	1
L-Leucine	Leu	¹³ C ₆ , 99%	5	1
L-Phenylalanine	Phe	ring- ¹³ C ₆ , 99%	100	1
Thymine	Т	1,3- ¹⁵ N ₂ , 98%	20	1
L-Tryptophan	Trp	¹³ C ₁₁ , 99%	100	1
L-Tyrosine	Tyr	ring- ¹³ C ₆ , 99%	100	1
Vitamin B ₃	nicotinamide	¹³ C ₆ , 99%	5	1
Citric acid	CA	1,5,6-carboxyl- ¹³ C ₃ , 99%	10	2
Fumaric acid	FA	¹³ C ₄ , 99%	100	2
Indole-3-acetic acid	IAA	phenyl-13C ₆ , 99%	5	2
α-Ketoglutaric acid, disodium salt	α-KG	1,2,3,4- ¹³ C ₄ , 99% (CP 97%)	100	2
Sodium palmitate	16:0	U- ¹³ C ₁₆ , 98%	10	2
Sodium pyruvate	Pyr	¹³ C ₃ , 99%	100	2


Companion unlabeled standard mixes and kits may be available; please inquire.

Chemical purity (CP) is 98% or greater, unless otherwise indicated.

Usage Specifications

Criteria	Recommendation		
Before reconstitution:			
Storage	ambient temperature		
Recommended retest	2 years from date of manufacture		
Upon reconstitution:			
Storage	4°C		
Recommended retest	4 weeks		

Example Results






Figure. Example XICs for a subset of labeled QReSS metabolites measured in buffer and human plasma by RPLC-MRM/MS (positive ESI, QTRAP 6500+). In the buffer analysis, a vial of a working stock mix was stored in on autosampler rack (maintained at 4°C) and processed eight days apart. Such measurements, when performed and metric tracked routinely, would constitute a system suitability test. In the human plasma analysis, an aliquot of the working stock mix was matrix spiked and subjected to a metabolomics workflow. The relative signals to its corresponding endogenous analyte (not shown for simplicity) is within an order of magnitude and can be used for relative quantitative applications, with absolute quantitation likely involving calibration curves. For additional information and application demonstrations see CIL application note #46.

Application Note

Percy, A.J.; Proos, R.; Demianova, Z.; Backiel, K.; Ubhi, B.K. 2021. Standardizing quantitative metabolomics analyses through the QReSS kit. (CIL application note #49)

The unique feature of our media analysis workflow is the use of the QReSS standard mix. We spike this into media samples before extraction, which helps normalize variabilities in metabolite extraction efficiencies and combat matrix effects. The QReSS mix consists of 18 isotopelabeled metabolites, which spans a multitude of metabolic classes and mimics the largely diverse chemical composition of cell culture media. This serves as an excellent choice for an internal standard mix and helps to improve the accuracy, precision, and robustness of our cell culture method.

– Hari Kosanam, PhD

Associate Principal Scientist, Vaccine Process Development and Commercialization

Merck & Co. (USA)

Metabolite Extract Kits

Expanded qualification remains a necessity in metabolomics today. Validated identifications are desired that can be used to springboard in-depth profiling studies and selective quantifications in preclinical screening studies. To aid method benchmarking (e.g., through a credentialing approach) and improved qualification/quantification studies, CIL offers two types of metabolite cell extracts: yeast (*Pichia pastoris*, strain CBS 7435) and *E. coli* (K12 strain MG1655), as outlined in the overview table below. Upon careful and precise solubilization (see reconstitution guidelines below), 100s of metabolites are potentially observable in their U-13C and unlabeled form. The metabolites span a broad class range (see example identification tables), having linkage to various biochemical pathways (e.g., citrate and glyoxylate cycle, amino acid and nucleotide metabolism) and cellular/molecular processes (e.g., immune system, blood coagulation, DNA metabolism). Please refer to the sample results and references below for examples of the product-specific applications.

Overview

Catalog No.	Description	No. of Metabolites	Kit Contents
ISO1	Metabolite Yeast Extract Kits (U-13C, 98%)	100s	1 vial of U-13C metabolite yeast extract document package (user manual, CoA, SDS, product flyer)
MSK-CRED-DD-KIT	Credentialed <i>E. coli</i> Cell Extract Kit (dried down)	100s	 1 vial (blue cap) of ¹³C-labeled <i>E. coli</i> cell extract (dried down) 1 vial (yellow cap) of unlabeled <i>E. coli</i> cell extract (dried down) document package (user manual, CoA, SDS, product flyer)
MSK-CRED-KIT	Credentialed <i>E. coli</i> Cell Extract Kit (solution)	100s	 1 vial (black) of ¹³C-labeled <i>E. coli</i> cell extract (solution) 1 vial (red) of unlabeled <i>E. coli</i> cell extract (solution) document package (user manual, CoA, SDS, product flyer)





Click on the thumbnails or visit isotope.com/applications/ for more information.

Metabolite Yeast Extract Kit

Compositions

Tabulated are the routinely identified metabolites in the yeast extracts (U-¹³C, ISO1; unlabeled, ISO1-UNL). **Note**: Additional metabolites have been measured with alternate methods and analysis techniques. This includes cofactors (e.g., NMN, ADPR, NADPH) and coenzymes (e.g., acetyl and malonyl coenzyme A).

Amino Acids and Derivatives (L enantiomer where applicable)				
S-Adenosyl-homocysteine (SAH)	Glutamate (Glu)	Methionine (Met)		
Alanine (Ala)	Glutamine (Gln)	(±)-3-Methyl-2-oxovalerate (K-IVal)		
α-Aminoadipate (AAD)	Glycine (Gly)	Ornithine (Orn)		
Arginine (Arg)	Guanidineacetate (GAA)	Phenylalanine (Phe)		
Argininosuccinate (ASA)	Histidine (His)	Proline (Pro)		
Asparagine (Asn)	Homoserine (Hse)	Sarcosine (Sar)		
Aspartate (Asp)	Isoleucine (Ile)	Serine (Ser)		
Betaine (BET)	α-Ketoisovalerate (KIV)	Threonine (Thr)		
Citrulline (Cit)	Kynurenine (KYN)	Tryptophan (Trp)		
Cystathionine (CYS)	Leucine (Leu)	Tyrosine (Tyr)		
Dihydroxyisovalerate (DIHV)	Lysine (Lys)	Valine (Val)		
Nucleobases, Nucleosides, and Nucleotides				
Adenine (Ade)	Cytidine triphosphate (CTP)	Inosine (Isin)		
Adenosine (Asin)	Deoxyadenosine monophosphate (dAMP)	Inosine monophosphate (IMP)		
Adenosine diphosphate (ADP)	5'-Deoxy-5'-methylthioadenosine (MTAP)	5-Methyluridine (m ⁵ U)		
Adenosine monophosphate (AMP)	Guanine (Gnin)	Pseudouridine (PsU)		
Adenosine triphosphate (ATP)	Guanosine (Gsin)	Uridine (Uri)		
Cyclic adenosine monophosphate (cAMP)	Guanosine diphosphate (GDP)	Uridine diphosphate (UDP)		
Cyclic guanosine monophosphate (cGMP)	Guanosine monophosphate (GMP)	Uridine monophosphate (UMP)		
Cytidine monophosphate (CMP)	Guanosine triphosphate (GTP)	Uridine triphosphate (UTP)		
Organic Acids				

Companion unlabeled standard mixes and kits may be available; please inquire.

Chemical purity (CP) is 98% or greater, unless otherwise indicated.

cis-Aconitate (cis-Ac)	DL-2-Hydroxyglutarate (2-HG)	Malate (Mal)
Citrate (CA)	Isocitrate (Iso)	Pyruvate (Pyr)
Fumarate (Fum)	α-Ketoglutarate (AKG)	Succinate (Suc)
Gluconate (GA)	Lactate (Lac)	
Sugar and Sugar Phosphates (D enantion	ner where applicable)	
Dihydroxyacetone phosphate (DHAP)	Glucose (Glc)	6-Phosphogluconate (6PGA)
Erythritol (Erthrit)	Glucose-6-phosphate (G6P)	Ribose (RIB)
Fructose (Fuc)	Mannitol (Man-Ol)	Ribose-5-phosphate (R5P)
Fructose-1,6-bisphosphate (FBP)	Mannose (Man)	Sedoheptulose-7-phosphate (S7P)
Fructose-6-phosphate (F6P)	Mannose-6-phosphate (M6P)	Trehalose (TRE)
Galactose (Gal)	2-Phosphoglycerate (2PG)	
Vitamins and Coenzymes		
Biotin (B ₇)*	Nicotinamide (NAM)	Nicotinamide adenine dinucleotide, reduced (NADH)
Choline (CHOL)	Nicotinamide adenine dinucleotide, oxidized (NAD+)	Nicotinamide adenine dinucleotide phosphate, oxidized (NADP+)
Other Small Molecules		
Glutamylcysteine (Glu-Cys)	Glutathione, reduced (GSH)	
Glutathione, oxidized (GSSG)	Mevalonate (MVA)	

^{*}Identified in ISO1-UNL only.

Usage Specifications

Extract Type	ISO1	ISO1-UNL	
No. of uses	~50*		
Before reconstitution:			
Storage	-80°C; protect from light		
Recommended retest	every 12 months		
Expiration	4 years from date of manufacture		
Upon reconstitution:			
Storage	4°C		
Recommended retest	4 weeks		

^{*}Conservative estimate, with numbers varying depending on method and application.

Solution Preparation Procedure

- 1. Reconstitute the yeast extract (ISO1 or ISO1-UNL) with 2 mL water, 2 mL 50% methanol, or 2 mL isopropanol.
- 2. Vigorously shake by hand with intermittent high-speed vortexing (2 minute minimum).
- 3. Centrifuge at 20°C for 5 min at 4000 rcf.
- 4. The clear standard solution can then be diluted (1/10 v/v) for direct use or prepared further for calibration and matrix addition.



A video demonstration of this procedure can be obtained here, as well as in the resources section of the Metabolite Yeast Extracts application page (see isotope.com/applications/metabolic research/metabolomics mixtures and kits).




Figure. Representative XICs of 13 C-labeled cofactors measured by HILIC-MS (positive ESI, Q Exactive HF). Shown in **A** are the routinely identified metabolites and **B** the additionally observed cofactors. Procedurally, the extract was reconstituted in 2 mL water then diluted 1:10 before analysis. The HILIC separation utilized an Acquity UPLC BEH Amide column ($100 \times 2.1 \text{ mm}$, $1.7 \mu m$ particles).

Example References

Mairinger, T.; Weiner, T.; Hann, S.; et al. **2020**. Selective and accurate quantification of *N*-acetylglucosamine in biotechnological cell samples via GC-MS/MS and GC-TOFMS. *Anal Chem, 92(7), 4875-4883*.

Rusz, M.; Rampler, E.; Keppler, B.K.; et al. **2019**. Single spheroid metabolomics: optimizing sample preparation of three-dimensional multicellular tumor spheroids. *Metabolites*, *9*(*12*), 304.

Galvez, L.; Rusz, M.; Schwaiger-Haber, M.; et al. **2019**. Preclinical studies on metal based anticancer drugs as enabled by integrated metallomics and metabolomics. *Metallomics*, *11(10)*, 1716-1728.

Sullivan, M.R; Danai, L.V.; Lewis, C.A.; et al. **2019**. Quantification of microenvironmental metabolites in murine cancers reveals determinants of tumor nutrient availability. *Elife*, *8*, e44235.

Demarest, T.G.; Truong, G.T.D.; Lovett, J.; et al. **2019**. Assessment of NAD⁺ metabolism in human cell cultures, erythrocytes, cerebrospinal fluid and primate skeletal muscle. *Anal Biochem, 572,* 1-8.

Hermann, G.; Schwaiger, M.; Volejnik, P.; et al. 2018. ¹³C-labelled yeast as internal standard for LC-MS/MS and LC high resolution MS-based amino acid quantification in human plasma. *J Pharm Biomed Anal, 155,* 329-334.

Guijas, C.; Montenegro-Burke, J.R.; Domingo-Almenara, X.; et al. **2018**. METLIN: A technology platform for identifying knowns and unknowns. *Anal Chem, 90(5),* 3156-3164. Si-Hung, L.; Causon, T.J.; Hann, S. **2017**. Comparison of fully wettable RPLC stationary phases for LC-MS-based cellular metabolomics. *Electrophoresis, 38(18),* 2287-2295.

Credentialed E. coli Extract Kit

Compositions

Condensed table of example metabolites detected in the E. coli extract samples (see PMID: 29256075 for application details). The spread of logPs (source: ALOGPS) highlights the breadth of physicochemical diversity.

Compound (Abbrev. or Alt. Name)	LogP	Metabolite Class
Adenine (Ade)	-0.38	Nucleobase
Adenosine monophosphate (AMP)	-2.30	Nucleotide
Biotin (vitamin B ₇)	0.30	Vitamin
Coenzyme Q10 (CoQ10; ubiquinone)	9.94	Coenzyme
Diacylglycerol (DG) 17:0/17:0	10.16	Lipid
Diacylglycerol (DG) 18:0/18:2	10.28	Lipid
Diacylglyercol (DG) 18:1/18:1	10.26	Lipid
Elaidic acid (EA; trans 18:1)	7.68	Fatty acid
Glutathione (GSH)	-2.70	Peptide
Oleic acid (OA; 18:1)	7.68	Fatty acid
Palmitic acid (PAL; 16:0)	7.23	Fatty acid
Palmitoyl CoA (PAL-CoA)	2.35	Coenzyme
Phosphatidylcholine (PC) 18:2/18:2	5.68	Lipid
Phosphatidylethanolamine (PE) 16:1/16:1	7.89	Lipid
Phosphatidylethanolamine (PE) 18:1/18:1	8.81	Lipid
Phenethylamine (PEA)	1.41	Neurotransmitter
Proline (Pro)	-2.70	Amino acid
Phosphatidylserine (PS) 18:1/18:1	4.96	Lipid
Retinol (vitamin A)	6.38	Vitamin
Stearic acid (STE; 18:0)	8.02	Fatty acid
Thymidine (Thd)	-1.30	Nucleoside
Tryptophan (Trp)	1.04	Amino acid
Uracil (U)	-1.20	Nucleobase
Uridine (Uri)	-1.80	Nucleoside

Usage Specifications

Criteria	Recommendation
No of sample runs	≤20 (with 5 µL LC-MS injections)
Storage	-80°C; protect from light, air, and moisture
Recommended retest	2 years

Solution Preparation Procedure

- 1. Reconstitute the dried-down extracts in 100 µL of ACN:water (e.g., 1:1 v/v). Note: Alternate solvent ratios are possible without impact to the credentialing approach. Changing the ratio to increase ACN or water will simply increase the coverage of organic or watersoluble metabolites, respectively.
- 2. Sonicate briefly the solutions (equates to 40 µL/mg of E. coli extract) then centrifuge (e.g., for 15 min at 13,000 RPM and 4°C) and incubate overnight at 4°C. This results in a clear solution.
- 3. Mix aliquots of the unlabeled and ¹³C-labeled extracts into autosampler vials at distinct ratios (e.g., 1:2 and 1:1 v/v). Vortex briefly.
- 4. Load the vials onto an autosampler rack for the benchmarking optimizations.




Figure. Credentialing benchmark approach for untargeted metabolomics. **A** Schematic of the credentialing workflow (abridged from **PMID: 25160088**) used to discriminate features of biological origin from contaminants and artifacts. The filtering criteria is based on feature ratios, with I being the intensity. **B** Credentialed metabolite measured on a Q Exactive Plus. Shown are the LC-MS details and spectra for glutathione ($C_{10}H_{17}N_3O_6S$). The M+0 in the MS survey scan is at m/z 308.0914, while the M+U is at 318.1249.

Example References

Sindelar, M.; Patti, G.J. 2020. Chemical discovery in the era of metabolomics. J Am Chem Soc, 142(20), 9097-9105.

Wang, L.; Naser, F.J.; Spalding, J.L.; et al. 2019. A protocol to compare methods for untargeted metabolomics. Methods Mol Biol, 1862, 1-15.

Spalding, J.L.; Naser, F.J.; Mahieu, N.G.; et al. **2018**. Trace phosphate improves ZIC-pHILIC peak shape, sensitivity, and coverage for untargeted metabolomics. *J Proteome Res*, *17(10)*, 3537-3546.

Naser, F.J.; Mahieu, N.G.; Wang, L.; et al. **2018**. Two complementary reversed-phase separations for comprehensive coverage of the semipolar and nonpolar metabolome. *Anal Bioanal Chem, 410(4),* 1287-1297.

Mahieu, N.G.; Patti, G.J. 2017. Systems-level annotation of a metabolomics data set reduces 25,000 features to fewer than 1,000 unique metabolites. *Anal Chem, 89(19),* 10397-10406.

Benton, H.P.; Ivanisevic, J.; Mahieu, N.G.; et al. 2015. Autonomous metabolomics for rapid metabolite identification in global profiling. Anal Chem, 87(2), 884-891.

Mahieu, N.G.; Huang, X.; Chen, Y.; et al. 2014. Credentialing features: a platform to benchmark and optimize untargeted metabolomic methods. Anal Chem, 86(19), 9583-9589.

Companion unlabeled standard mixes and kits may be available; please inquire.

Chemical purity (CP) is 98% or greater, unless otherwise indicated.

IROA® Biochemical Quantitation Kits

IROA's quantitative assay kits can be used to study physiological stressors (such as disease or environmental factors), phenotypically distinct cell lines, biomarkers, systems biology, and flux in a wide variety of cell populations and biological samples (see references below for background and application examples). These kits utilize 13C-labeled energy sources at reduced isotopic enrichment (e.g., 95% U-13C and 5% U-13C D-glucose) in their respective control and experimental populations, which helps with the predictable and distinguishable detection of MS-based isotopic distributions. These distributions can then be used to: (i) differentiate biological signals from artifacts, (ii) calculate accurate molecular formulae, and (iii) determine relative concentrations of the metabolites of biological origin. The details of the metabolic profiling kits (bacterial and mammalian) available off-the-shelf are outlined below. Also specified below is the kit usage specifications, a protocol schematic, and example references.



isotope.com/applications/ for

Overview

Catalog No.	Description	Protocol Type	Kit Contents
IROA-200-50	IROA 200 Kit for Bacterial Metabolic Profiling	Basic	 Control medium: M9 minimal medium containing 95% U-¹³C D-glucose Experimental medium: M9 minimal media containing 5% U-¹³C D-glucose Control component mix: 95% U-¹³C amino acids Experimental component mix: 5% U-¹³C amino acids Guides and tools: user manual, ClusterFinder™ software, statistical analysis package
IROA-300-250	IROA 300 Kit for Mammalian Metabolic Profiling	Basic, Phenotypic, or Fluxomic	IROA PHENO-95-300 (for control cell population labeling, see contents below) IROA FLUX-05-300 (for experimental cell population labeling, see contents below) Note: The kits can be used independently (for phenotypic or fluxomic profiling) or combined (for basic profiling).
IROA-PHENO-95-300	IROA 300 Kit for Phenotypic Metabolic Profiling	Phenotypic	 Control medium: 95% U-¹³C D-glucose Control component mix: 95% U-¹³C amino acid mix and 95% U-¹³C yeast extract Vitamins: EBSS/RPMI 1640 Guides and tools: user manual, ClusterFinder software, statistical analysis package
IROA-FLUX-05-300	IROA 300 Kit for Fluxomic Metabolic Profiling	Fluxomic	 Experimental medium: 5% U-¹³C D-glucose Experimental component mixes: 5% U-¹³C amino acid mix and 5% U-¹³C yeast extract Vitamins: EBSS/RPMI 1640 Guides and tools: user manual, ClusterFinder software, statistical analysis package

IROA is a registered trademark of IROA Technologies.

Usage Specifications

Kit Materials	Storage	Stability
Liquid media	4°C; protect from light	2 years from date of manufacture
Neat component mixes	-20°C; in darkness	2 years from date of manufacture
Reconstituted component mixes	-20°C; in darkness	1 year

The sample numbers per kit are predicated on the method. Noted here are the minimal number of sample uses based on the method supplied in the user manual. Procedurally for the IROA 300s, this involves a six-well plate growth, 2.5 generations per passage, and 3 mL of media/well, with five cell doublings implemented to ensure full label incorporation.

Kit Type	Implemented Protocol	Minimal No. of Uses
IROA-200-50	Basic	48 (control and experimental samples)
IROA-300-250	Basic	72 (control and experimental samples)
IROA-PHENO-95-300	Phenotypic	72 (control samples)
IROA-FLUX-05-300	Fluxomic	72 (experimental samples)

Schematic of Protocols




Figure. MS spectra for L-arginine $(C_6H_{14}N_4O_2)$ measured by LC-MS under three types of IROA protocols. Briefly, the cell populations are grown with isotopically labeled carbon sources, which includes D-glucose for: (1) control and experimental samples in the basic protocol, (2) control samples only in the phenotypic protocol, and (3) experimental samples only in the fluxomic protocol (tracers added after harvest at 95 or 99% 13 C). **Note:** The control or tracer signals (at 95% 13 C in basic and phenotypic protocols or 99% 13 C in flux) are illustrated in blue and the experimental signals (at 5% 13 C in basic and flux protocols or natural abundance in phenotypic) in green or black.

Example References

Carey, J; Nguyen, T; Korchak, J.; et al. **2019**. An isotopic ratio outlier analysis approach for global metabolomics of biosynthetically talented actinomycetes. *Metabolites*, *9*(*9*), 181. Vlahakis, C.; Hazebroek, J.; Beecher, C.; et al. **2019**. Isotopic ratio outlier analysis improves metabolomics prediction of nitrogen treatment in maize. *Phytochemistry*, *164*, 130-135. Rahman, M.A.; Babar, A; Akond, M.; et al. **2017**. LC-HRMS Based non-targeted metabolomic profiling of wheat (*Triticum aestivum l*.) under post-anthesis drought stress. *Am J Plant Sci*, *8*(*12*), 3024-3061.

Viant, M.R.; Kurland, I.J.; Jones, M.R., et al. 2017. How close are we to complete annotation of metabolomes? Curr Opin Chem Biol, 36, 64-69.

Clendinen, C.S; Stupp, G.S.; Ajredini, R.; et al. **2015**. An overview of methods using (13)C for improved compound identification in metabolomics and natural products. *Front Plant Sci.* 6. 611.

Edison, A.S.; Clendinen, C.S.; Ajredini, R.; et al. **2015**. Metabolomics and natural-products strategies to study chemical ecology in nematodes. *Integr Comp Biol, 55(3),* 478-485. de Jong, F.A.; Beecher, C. **2012**. Addressing the current bottlenecks of metabolomics: Isotopic Ratio Outlier AnalysisTM, an isotopic-labeling technique for accurate biochemical profiling. *Bioanalysis, 4(18),* 2303-2314.

PeptiQuant™ Plus Assay Kits

for each kit type are listed in the overview below.

Researchers in academia and life science industries continue to implement a bottom-up MS-based workflow for protein biomarker screening. Biomarker verification/validation requires absolute quantification of surrogate peptides in a sample matrix, a requirement that is best achieved using well-characterized stable isotope-labeled standards. To ensure robust quantitative measurements, QC checks should be routinely performed. CIL offers a collection of PeptiQuant Assay Kits (from MRM Proteomics Inc.) for QC and biomarker assessment using bottom-up LC-MS/MS methodologies. The QC kits are designed to evaluate the performance of an LC-MS platform, either alone or in combination with a human plasma proteomics workflow (see corresponding flyer for panel details). The biomarker assessment kits (BAKs) are intended to help researchers screen target panels of candidate protein disease biomarkers in human or mouse plasma samples (see corresponding flyer for panel details). The current platform-specific offerings

MRM Proteomics mrmproteomics.com

Overview

QC Kits					
Catalog No.	Description	Kit Contents	No. of Peptides	Unit Size	Optimized Instrument
LCMSP-QC	PeptiQuant Plus Human Plasma Daily QC Kit	13C/15N-labeled peptide mix USB (e.g., user manual, acquisition and analysis files)	35	10, 20, or 50 injections	 6490/6495 QqQ QTRAP® 6500 Q Exactive™ Plus
WFPK	PeptiQuant Plus Human Plasma Workflow QC Kit	13C/15N-labeled peptide mix unlabeled peptide mix trypsin BSA human plasma USB (e.g., user manual, acquisition and analysis files)	35	1 or 2 runs	• 6490/6495 QqQ • QTRAP 6500 • Q Exactive Plus

Biomarker As	Biomarker Assessment Kits (BAKs)				
Catalog No.	Description	Kit Contents	No. of Peptides	Unit Size	Optimized Instrument
BAK-125	PeptiQuant Plus Human Plasma Proteomics Kit	13C/15N-labeled peptide mix unlabeled peptide mix trypsin BSA USB (e.g., user manual, acquisition and analysis files)	125	20, 50, or 100 samples	• 6490/6495 QqQ • QTRAP 6500 • Q Exactive Plus • Xevo TQ-XS
M-BAK-125*	PeptiQuant Plus Mouse Plasma Proteomics Kit	13C/15N-labeled peptide mix unlabeled peptide mix trypsin BSA USB (e.g., user manual, acquisition and analysis files)	125	20, 50, or 100 samples	• 6490/6495 QqQ • QTRAP 6500 • Q Exactive Plus • 6545 Q-TOF
BAK-270	Expanded PeptiQuant Plus Human Plasma Proteomics Kit	13C/15N-labeled peptide mix unlabeled peptide mix trypsin BSA USB (e.g., user manual, acquisition and analysis files)	270	100 samples	• 6490/6495 QqQ • QTRAP 6500 • Q Exactive Plus





Click on the thumbnails or visit isotope.com/applications/ for more information

PeptiQuant is a trademark of MRM Proteomics Inc.

Usage Specifications

3 1	
Before reconstitution:	
Storage	-80°C
Recommended retest	6 months from date of manufacture

Note: The dried-down trypsin and BSA materials supplied with the BAKs are not stable at room temperature and must be stored as above. Trypsin is to be prepared immediately prior to use and stored on ice until dispensed.

^{*}Alternate sets of 125 target proteins are available (see PeptiQuant Plus BAK flyer for details).




Figure. Bottom-up LC-MRM/MS data for two representative, low abundance targets from a PeptiQuant Plus BAK-270 analysis. Illustrated are Skyline screenshots for peptide DVIVHPLPLK (from Adhesion G protein-coupled receptor F5, gene ADGRF5, UniProtKB Q8IZF2) in **A** and peptide LLEYDTVTR (from Adipocyte plasma membrane-associated protein, gene APMAP, UniProtKB Q9HDC9) in **B**. The tabs refer to the monitored panel, two example XICs (for curve calibrant level F and a pooled human plasma sample analysis), the standard curve (with calibrant levels marked with gray squares, QC samples with green diamonds, experimental sample with red arrow), and the results grid. The protein concentration in the pooled human plasma sample was determined to be 321.2 ng/mL for ADGRF5 and 439.9 ng/mL for APMAP. Results courtesy of MRM Proteomics.

Example References

Percy, A.J.; Borchers, C.H. 2021. Detailed method for performing the ExSTA approach in quantitative bottom-up plasma proteomics. Methods Mol Biol, doi: 10.1007/978-1-0716-1024-4.

Ayton, S.; Janelidze, S.; Roberts, B.; et al. 2021. Acute phase markers in CSF reveal inflammatory changes in Alzheimer's disease that intersect with pathology, APOE &4, sex and age. Prog Neurobiol, 101904.

Gaither, C.; Popp, R.; Mohammed, Y.; et al. 2020. Determination of the concentration range for 267 proteins from 21 lots of commercial human plasma using highly multiplexed multiple reaction monitoring mass spectrometry. Analyst, 145(10), 3634-3644.

Michaud, S.A.; Sinclair, N.J.; Petrošová, H.; et al. 2018. Molecular phenotyping of laboratory mouse strains using 500 multiple reaction monitoring mass spectrometry plasma assays

Orti, V.; Mertens, B.; Vialaret, J.; et al. 2018. Data from a targeted proteomics approach to discover biomarkers in saliva for the clinical diagnosis of periodontitis. Data Brief, 18, 294-299.

Mohammed, Y.; Pan, J.; Zhang, S.; et al. 2018. ExSTA: External standard addition method for accurate high-throughput quantitation in targeted proteomics experiments. Proteomics Clin Appl, 12(2), 1600180.

The PeptiQuant Plus Platform Performance Kit has proven to be a vital component of our everyday quality assurance that enables us to deliver high-quality targeted proteomics data in an accurate and timely manner. This kit has a 'dilute and shoot' operation and comes with vendorspecific LC-MRM/MS parameters and a Skyline analysis file for quick input and results output. Altogether, the performance kit is an excellent means to rapidly assess LC-MS performance that should become a routine staple in a proteomic user's toolbox.

Associate Principal Scientist at AstraZeneca Gothenburg (Sweden)

FeptiQuant Plus Assay Kits contain all the essential materials, including the standards and methods, for performing absolute protein quantification by LC-MRM/MS in a standardized way. The standard protocol helped us reduce the assay development time, while improve the reproducibility and precision of multiplex protein quantification. In addition to the biomarker assessment kits, the quality control kits enable the instrument performance and assay reproducibility to be monitored and assessed, which ultimately provided us confidence in the reliability of the quantification results.

- Elaine Wong, PhD

Scientific Officer at Queen Mary Hospital, Fu Lam (Hong Kong)

INLIGHT® Glycan Tagging Kit

Glycans participate in a large number of cellular, molecular, and biological processes and are implicated in a number of diseases (e.g., Alzheimer's, cancer). To aid the identification and relative quantification of glycans by LC-MS, CIL offers the innovative INLIGHT® (Individuality Normalization when Labeling with Isotopic Glycan Hydrazide Tag) glycan-tagging kit. This kit employs both natural (NAT – $^{12}C_6$) and stable isotope-labeled (SIL – $^{13}C_6$) phenyl 2-GPN reagents in the hydrazide derivatization of free *N*-glycans. While this tagging strategy was developed for *N*-glycans, it has also been adapted to sample analysis of *O*-glycan and heparin oligomer profiles. The INLIGHT kit contains five vials of NAT reagent and five vials of SIL reagent, which in total provides sufficient tagging for approximately 125 relative quantification experiments. The user manual provides step-by-step instructions for executing the modified INLIGHT strategy using maltoheptaose and fetuin A as examples. Data processing and analysis of derivatized glycans can be facilitated in GlycoHunter or Skyline.



Click on the thumbnails or visit isotope.com/applications/ for more information

Overview

Catalog No.	Description	Kit Contents
GTK-1000	INLIGHT® Glycan Tagging Kit	• Light phenyl 2-GPN INLIGHT reagent (5 × 0.25 mg)
		• Heavy phenyl 2-GPN INLIGHT reagent (${}^{13}C_6$; 5 × 0.25 mg)
		• Unlabeled maltoheptaose (5 × 10 μg)

Usage Specifications

Criteria	Recommendation
No. of Uses	25 per vial
Before reconstitution:	
Storage	ambient temperature; protect from light and moisture
Recommended retest	5 years from date of manufacture




Figure. LC-MS examples of INLIGHT derivatized maltoheptaose and *N*-linked glycans from bovine fetuin.

- **A**. Extracted ion chromatogram of [NAT + H]⁺ maltoheptaose.
- **B**. Mass spectrum at 10 minutes illustrating NAT and SIL derivatized maltoheptaose, maltohexaose, maltopentaose, and maltotetraose.
- **C**. Extracted ion chromatogram of a fetuin-based glycan at m/z 1012.3980 corresponding to [NAT + 2H + (Fuc)₁(Gal)₂(GlcNAc)₄(Man)₃]²⁺.
- ${\bf D}$. Average mass spectrum of a fetuin analysis showing the NAT and SIL derivatized N-linked glycan pairs in the LC window 8-11 min.

Example References

Kalmar, J.G.; Garrard, K.P.; Muddiman, D.C. **2021**. GlycoHunter: An open-source software for the detection and relative quantification of INLIGHT®-labeled *N*-linked glycans. *J Proteome Res, 20(4)*, 1855-1863.

Kalmar, J.G.; Butler, K.E; Baker, E.S.; et al. **2020**. Enhanced protocol for quantitative *N*-linked glycomics analysis using Individuality Normalization when Labeling with Isotopic Glycan Hydrazide Tags (INLIGHT).* *Anal Bioanal Chem, 412(27), 7569-7579*.

King, S.R.; Hecht, E.S.; Muddiman, D.C. 2018. Demonstration of hydrazide tagging for O-glycans and a central composite design of experiments optimization using the INLIGHT® reagent. Anal Bioanal Chem, 410(5), 1409-1415.

Loziuk, P.L.; Hecht, E.S.; Muddiman, D.C. **2017**. *N*-linked glycosite profiling and use of Skyline as a platform for characterization and relative quantification of glycans in differentiating xylem of *Populus trichocarpa*. *Anal Bioanal Chem*, *409(2)*, 487-497.

Hecht, E.S.; Scholl, E.H.; Walker, S.H.; et al. **2015**. Relative quantification and higher-order modeling of the plasma glycan cancer burden ratio in ovarian cancer case-control samples. *J Proteome Res*, *14*(10), 4394-4401.

Chemical purity (CP) is 98% or greater, unless otherwise indicated.

O & As

Listed below are a series of general Q&As for CIL's mixes and kits. Product-specific FAQs can be located at isotope.com/applications/ under their corresponding product application page (e.g., Metabolic Research → Metabolomics Mixtures and Kits → QReSS Kits).

What are the advantages of selecting an off-the-shelf mix vs. a collection of individual isotope standards for self-mixing?

The MSK/NSK products are formulated to exacting standards following detailed batch records developed from over 20 years of formulation experience. Following production, randomly selected vials are analyzed to ensure both accuracy and consistency. The astute attention and process control generates exceptional vial-tovial and lot-to-lot reproducibility. This high reproducibility return is one of the merits of utilizing CIL prepared mixes. Additional advantages are:

- reduced development time and cost;
- enhanced data quality;
- ease of user implementation; and
- improved confidence in analytical results.

The majority of the mixes are ¹³C- and/or ¹⁵N-labeled. What are the advantages of ¹³C/¹⁵N-labeling vs. D-labeled in MS measurements?

The nature of the stable isotope can potentially impact the preanalytical (e.g., storage and handling) and analytical (e.g., sample preparation and processing) phases of an experiment. In comparison to D labels, ¹³C/¹⁵N labels can have:

- improved isotope stability;
- negligible isotope scrambling issues;
- conserved chromatographic elution (relative to its unlabeled standard); and
- · heightened analytical reliability.

In cases where D-labeled compounds were selected, the D-label is located at nonexchangeable positions and was stability tested for preservation as well as product application. For further background, please refer to our technical note (on p. 9 of our Stable Isotope Standards for Clinical Mass Spectrometry catalog) that describes the benefits of ¹³C vs. D standards in MS-based studies.

How often should QC measurements be performed, and what mixes are most suitable?

Prior to first use. OCs should be conducted several times in succession prior to sample analysis to establish baseline performance and intervention limits. Once established, QCs should be performed routinely (i.e., before, during, and after sample analysis) to monitor the effectiveness of the analytical method and instrument platform over time. While all of CIL's mixes could be used in these types of operations, the ones ideally suited for these type of measurements are the metabolomic QC kits (see MSK-QC-KIT and MSK-QReSS-KIT) and the proteomic QC kits (see PeptiQuant Plus daily and workflow QC kits).

In a targeted quantitative application, how many MRM transitions should be monitored per analyte?

Ideally, a minimum of two MRM transitions per compound (be it metabolite or peptide) should be targeted. This will allow for a quantifier and qualifier(s) assignment, with the ion ratios serving as an additional metric for performance qualification. Nonetheless, this may not be possible for all analytes given the compound's fragmentation chemistry, the employed MS/MS parameters, and the mode of operation utilized. Empirical experiments with isotope-labeled standards should therefore be first performed to optimize the MS/MS parameters prior to conducting the quantitative study with precious experimental samples.

Can custom mixes or add-on vials be formulated?

Yes, we have the ability to customize. We would first review feasibility and then provide a quotation on your specific mix. To start this process, please provide the necessary details on this custom mix request form or contact your local sales representative.

Please visit isotope.com for a complete list of isotope-labeled compounds.

Research products are distributed and sold worldwide via our extensive network.

CIL's distributor listing is available at isotope.com.

To request a quotation or place an order:

North America: 1.978.749.8000 | 1.800.322.1174 | cilsales@isotope.com International: +1.978.749.8000 | intlsales@isotope.com Fax: 1.978.749.2768 | isotope.com