



CIL

Cambridge Isotope Laboratories, Inc.
isotope.com

Eurisotop

A CAMBRIDGE ISOTOPE LABORATORIES COMPANY

RESEARCH PRODUCTS

Standards for qNMR

qNMR Standards for External and Internal Referencing



Quantitative ¹H-NMR (qNMR) continues to be utilized with much success in the pharmaceutical, chemical and food industries and in many facets of academic research. Regardless of the application, all qNMR methods require a calibration signal whose integrated signal intensity originates or is traceable to a known number of protons. Calibration for qNMR is made using either internal or external referencing methods. External methods rely on the use of a standard solution packaged in a defined NMR tube or capillary to obtain an integral that can be used for sample quantification, whereas internal methods rely on the use of a known amount of standard that is co-dissolved in the sample itself.

External Calibration Standards

CIL is pleased to offer external calibration standards for qNMR. The standards are formulated using CIL's high-quality DMSO-d₆ and benzoic acid from NIST (SRM 350(b)), a standard reference material for acidometry. Both 5 mM and 15 mM benzoic acid concentrations are available. The concentration and associated expanded uncertainty of the benzoic acid has been accurately determined using metrological techniques and verified using qNMR. The ¹H-NMR spectrum of benzoic acid in DMSO-d₆ is presented in Figure 1.


CIL is currently offering these standards in presealed NMR tubes. Please see the information below for details regarding NMR tubes and fill volumes. Other NMR tubes and concentrations may be available upon request.

qNMR Standard for External Referencing

| Catalog No. | Description* | NMR Tube** | Part No. | Fill Volume |
|-------------|---|-------------|--------------------------|-------------|
| DLM-9491A | 5 mM Benzoic acid in DMSO-d ₆ | 1.7 mm O.D. | Bruker Part No. Z106462 | 50 µL |
| DLM-9491B | 5 mM Benzoic acid in DMSO-d ₆ | 3 mm O.D. | Wilmad Part No. 335-PP-9 | 160 µL |
| DLM-9491C | 5 mM Benzoic acid in DMSO-d ₆ | 5 mm O.D. | Wilmad Part No. 528-PP-8 | 750 µL |
| DLM-7061A | 15 mM Benzoic acid in DMSO-d ₆ | 1.7 mm O.D. | Bruker Part No. Z106462 | 50 µL |
| DLM-7061B | 15 mM Benzoic acid in DMSO-d ₆ | 3 mm O.D. | Wilmad Part No. 335-PP-9 | 160 µL |
| DLM-7061C | 15 mM Benzoic acid in DMSO-d ₆ | 5 mm O.D. | Wilmad Part No. 528-PP-8 | 750 µL |

* The benzoic acid concentration and associated uncertainty are reported.

** All tubes are flame-sealed to ensure longevity.



Other NMR fill volumes and tubes are available. Please inquire.



Eurisotop, Parc des Algorithmes, route de l'orme, 91190 Saint Aubin | France

tel: +33 1 69 41 97 98

fax: +33 1 69 41 93 52

+49 (0) 681 99 63 338 (Germany)

www.eurisotop.com

Internal Calibration Standards

The internal reference method commonly gives errors of <1% and is considered to be the most accurate and reproducible method available to obtain quantitative $^1\text{H-NMR}$ spectra. Unfortunately, the reference standard is typically weighed into each sample solution, an action that requires time and effort, and has been reported to be the largest source of error with this method.

CIL is pleased to offer a ready-to-use DMSO-d_6 solution containing a known amount of benzoic acid for internal referencing. Because this solution is preformulated, the user does not need to weigh a

standard material. The elimination of this step will reduce effort and time in sample preparation and also may bring about more accurate results than if the user performs this formulation. To use this product, the sample must be soluble in DMSO-d_6 , physically and chemically inert toward benzoic acid and stable in acidic pH. Ideally, there will be no resonances from the sample in the region of benzoic acid aromatic protons (7.4-8.1 ppm), HOD (~3 ppm but is variable) and DMSO-d_5 (2.5 ppm). The benzoic acid concentration with associated uncertainty is presented on the certificate of analysis.


qNMR Standard for Internal Referencing

| Catalog No. | Description | Ampoule | Comments |
|-------------|---|---------|--|
| DLM-9491D | 5 mM Benzoic acid in DMSO-d_6 | 1 g | The benzoic acid concentration and associated uncertainty is reported. |
| DLM-7061D | 15 mM Benzoic acid in DMSO-d_6 | 1 g | The benzoic acid concentration and associated uncertainty is reported. |

CIL Formulation Procedure

The procedure that CIL uses to formulate qNMR external calibration reference standard bulk solutions allows for the expanded uncertainty of the concentration of the calibration standard (e.g., benzoic acid) to be determined. Traceability to SI is maintained through the use of weight sets with

calibration traceable to NIST and laboratory balances with NIST-traceable calibration certificates, maintaining an unbroken chain of calibration to the kilogram. The factors contributing to the uncertainty of the benzoic acid concentration¹ is shown in below.



Cause-and-effect diagram of factors contributing to the uncertainty of the benzoic acid concentration in the qNMR standard formulation.

Reference

1. EURACHEM CITAC Guide CG 4, "Quantifying Uncertainty in Analytical Measurement," Third Edition, QUAM:2012

