



CIL

Cambridge Isotope Laboratories, Inc.  
isotope.com

Eur isotop

A CAMBRIDGE ISOTOPE LABORATORIES COMPANY

RESEARCH PRODUCTS

# Metabolite Yeast Extract

## For Quality Control and Internal Standardization in Metabolomics



Cambridge Isotope Laboratories, Inc. (CIL) is pleased to offer a dilute-and-shoot  $^{13}\text{C}$ -labeled and unlabeled metabolite yeast extract for use as an internal standard in MS-based quantitative or profiling studies of various sample types.<sup>1-9</sup>

These extracts have been highly characterized by a number of methodologies and are amenable to a variety of research uses after simple reconstitution. The components in the extracts span broad metabolic classes (e.g., amino and organic acids, sugar phosphates, coenzymes), biochemical pathways (e.g., citrate and glyoxylate cycle, amino acid and nucleotide metabolism), and cellular/molecular processes (e.g., immune system, blood coagulation, DNA metabolism).

### Overview

#### Uses

- Targeted or untargeted, MS-based analysis
- Method and instrument QC
- Metabolite quantitation
- Biomarker discovery and verification

#### Benefits

- Reduce measurement uncertainty
- Improve precision and accuracy
- Enhance identification confidence
- Decrease development time and cost

### Description

Catalog No.	Description	Price
ISO1	Metabolite Yeast Extract ( $\text{U-}^{13}\text{C}$ , 98%)	\$500
ISO1-UNL	Metabolite Yeast Extract (unlabeled) <b>NEW!</b>	\$300

Dry extract of  $>2 \times 10^9$  *Pichia pastoris* cells (~15 mg).  
Produced by ISOTopic Solutions (isotopic-solutions.com).



### Specifications

Criteria	Note
Kit contents	<ul style="list-style-type: none"> <li>• <math>\text{U-}^{13}\text{C}</math>-labeled or unlabeled yeast extract</li> <li>• detailed user manual (on USB)</li> </ul>
Shipping	Dried extract in 15 mL Falcon™ tube on dry ice
Storage	-80°C, protected from light
Isotopic enrichment	>98% (via LC-MS for characteristic metabolites)
Chemical purity	99% (via LC-UV for characteristic metabolites)

### References

- Hermann, G.; Schwaiger, M.; Volejnik, P.; et al. **2018**.  $^{13}\text{C}$ -labelled yeast as internal standard for LC-MS/MS and LC high resolution MS-based amino acid quantification in human plasma. *J Pharm Biomed Anal*, 155, 329-334.
- Guijas, C.; Montenegro-Burke, J.R.; Domingo-Almenara, X.; et al. **2018**. METLIN: A technology platform for identifying knowns and unknowns. *Anal Chem*, 90(5), 3156-3164.
- Si-Hung, L.; Causon, T.J.; Hann, S. **2017**. Comparison of fully wetttable RPLC stationary phases for LC-MS-based cellular metabolomics. *Electrophoresis*, 38(18), 2287-2295.
- Schwaiger, M.; Rampler, E.; Hermann, G.; et al. **2017**. Anion-exchange chromatography coupled to high-resolution mass spectrometry: A powerful tool for merging targeted and non-targeted metabolomics. *Anal Chem*, 89(14), 7667-7674.
- Ortmayr, K.; Hann, S.; Koellensperger, G. **2015**. Complementing reversed-phase selectivity with porous graphitized carbon to increase the metabolome coverage in an on-line two-dimensional LC-MS setup for metabolomics. *Analyst*, 140(10), 3465-3473.
- Neubauer, S.; Chu, D.B.; Marx, H.; et al. **2015**. LC-MS/MS-based analysis of coenzyme A and short-chain acyl-coenzyme A thioesters. *Anal Bioanal Chem*, 407(22), 6681-6688.
- Klavins, K.; Drexler, H.; Hann, S.; et al. **2014**. Quantitative metabolite profiling utilizing parallel column analysis for simultaneous reversed-phase and hydrophilic interaction liquid chromatography separations combined with tandem mass spectrometry. *Anal Chem*, 86(9), 4145-4150.
- Klavins, K.; Chu, D.B.; Hann, S.; et al. **2014**. Fully automated on-line two-dimensional liquid chromatography in combination with ESI MS/MS detection for quantification of sugar phosphates in yeast cell extracts. *Analyst*, 139(6), 1512-1520.
- Neubauer, S.; Haberhauer-Troyer, C.; Klavins, K.; et al. **2012**.  $\text{U-}^{13}\text{C}$  cell extract of *Pichia pastoris* – a powerful tool for evaluation of sample preparation in metabolomics. *J Sep Sci*, 35(22), 3091-3105.



Euriso-Top, Parc des Algorithmes, route de l'orme, 91190 Saint Aubin | France

tel: +33 1 69 41 97 98

fax: +33 1 69 41 93 52

+49 (0) 681 99 63 338 (Germany)

www.eur isotop.com