



CIL

Cambridge Isotope Laboratories, Inc.
isotope.com

Eurisotop

A CAMBRIDGE ISOTOPE LABORATORIES COMPANY

¹⁸O WATER


¹⁸O Water

The Image of Quality



*With our partner, **ABX advanced biochemical compounds**, we provide complete PET chemistry solutions*





Committed to the quality the medical community demands

Consistently exceeding your expectations:

^{18}O Water


- $\geq 97\%$ ^{18}O enrichment
- $\leq 3\mu\text{mho/cm}$ conductivity
- State-of-the-art quality-control lab
- Produced in the world's largest stable isotope-separation facility in Xenia, OH

Consistent superior quality

Controlled process

Dependable yield

Reliable product



Cambridge Isotope Laboratories

Leadership

- With our partner, **ABX advanced biochemical compounds**, we provide **complete PET chemistry solutions**.
- **Preferred supplier** to the world's leading ^{18}F FDG producers.
- **World leader** in the separation of ^{18}O and ^{13}C isotopes for medical applications and the production of stable isotope-labeled research and medical compounds.



Commitment to PET

- Supporting the PET community with an **^{18}O water production facility** located in the US, with an annual production capacity of approximately 450 kg and the ability to expand to meet increased market demand.
- **With ABX, providing innovation in radiochemistry labeling strategies.**
- **For more than 35 years**, CIL products have contributed to advancements in drug discovery, environmental analysis, genomics and proteomics, physics, and medical diagnostic research.
- **CIL offers superior service to meet your needs:**
 - exceptional quality
 - convenient ordering
 - competitive pricing
 - dependable delivery
 - flexible packaging
 - extensive sales distribution network





To purchase ^{18}O water contact us at:
18Owater@isotope.com
or +1.978.749.8000

Cambridge Isotope Laboratories, Inc.
3 Highwood Drive, Tewksbury, MA 01876 USA
www.18Owater.com

